
Pokyny	CZ
Bedienungsanleitung	DE
Brugsvejledning	DK
Instrucciones	ES
Käyttöohje	FI
Mode d'emploi	FR
Instructions	GB
Οδηγίες	GR
Használati utasítás	HU
Manuale d'Uso	IT
Gebruiksaanwijzing	NL
Bruksanvisning	NO
Instrukcja obsługi	PL
Instruções	PT
Bruksanvisning	SE
Uputstvo za upotrebu	YU

Zesilovač signálu Bellman Visit 868, BE1510

Úvod

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Bellman Visit 868 Repeater, BE1510

Introduction

Thank you for choosing products from Bellman & Symfon.

The Bellman Visit 868 System consists of a number of radio transmitters and receivers. The transmitters detect different events in the surrounding area and transmit a radio signal to the receivers. The receivers pick up this signal and provide indications using light, sound and/or vibration.

The transmitter determines what type of light, sound or vibration should be displayed so that the reason for the indication is evident.

Read through the entire user manual first and then start to install the system.

Refer to the illustration of the Bellman Visit 868 System on the inside of the cover.

Getting started

Unpacking, installing and testing the unit

1. Open the battery cover (7). Connect the power supply unit to the socket (6), fit one back-up battery, either a 6LR61 alkaline or 6F22 lithium type battery, and close the battery cover.
2. Press the Test Button (2). The Bellman Visit 868 Repeater blinks with the LED (1) to indicate that the unit is transmitting a radio signal. The Bellman Visit 868 receivers will indicate the Door Signal.

3. Mount the Bellman Visit 868 Repeater on the wall either using the self-adhesive Velcro tape or suspended on the wallmount bracket (8) on the screw supplied. If the Velcro tape is used, the wall on which the transmitter is to be positioned can be cleaned using the enclosed wet wipe. The Repeater can also be placed on a level surface, e.g. a table, but avoid placing it on the floor.

Function

General

The Bellman Visit 868 Repeater, BE1510, is a repeater in the Bellman Visit 868 System for indoor use which receives a radio signal from a transmitter in the Bellman Visit 868 System and transmits this radio signal to the receivers in the system. The aim is to increase the range and obtain reception in locations in which it is not otherwise possible by using Bellman Visit 868 Repeaters.

The Repeater also operates as a multi-function transmitter, either via a pushbutton or via an input for an Exterior External Trigger. The Test Button and Exterior External Trigger can be detected separately and transmit different Signal Patterns to the receivers in the Bellman Visit 868 System depending on what has activated the Bellman Visit 868 Repeater. There is a range of options for connecting the Bellman Visit 868 Repeater to various applications. We recommend using a maximum of one Repeater per system.

Radio key

On delivery all Bellman Visit 868 units are tuned to the same Radio Key. If you have a neighbour with a similar system, you can change to different Radio Keys so that you do not affect each other's systems. If you use the Radio Key Switch (5) to change the Radio Key on this transmitter, you must also change all other units in your Bellman Visit 868 System to the same Radio Key. Refer to the user manual for the relevant unit.

Indicators and Signals

System indicators

The LED (1) blinks, when the Bellman Visit 868 Repeater transmits or receives radio signals.

Power supply

When the Repeater is connected to the power supply unit, the LED is constantly green, apart from when it is receiving or transmitting radio signals when it blinks.

When the Repeater is activated, the LED (1) lights up green. This means that the back-up battery is in good condition.

If the LED (1) blinks yellow, this means that the back-up battery is flat and must be changed. Only use a 6LR61 (alkaline) or 6F22 (lithium) type battery.

Troubleshooting in brief

Problem	Solution
Nothing happens when the transmitter is activated with the Test Button (2).	<ul style="list-style-type: none"> • Check that the power supply unit is connected correctly. • If the power supply unit is not connected, change the back-up battery. Only use an alkaline 6LR61 or a lithium 6F22 type battery. Make sure that the power supply unit is connected correctly.

The LED (1) blinks yellow when the Repeater is activated.	<ul style="list-style-type: none"> • Change the back-up battery. Only use an alkaline 6LR61 or a lithium 6F22 type battery. Make sure that the power supply unit is connected correctly.
The LED (1) blinks green when the Repeater is activated, but the receivers are not signalling.	<ul style="list-style-type: none"> • Check the battery in the receiver. • Check that the receivers are not placed too far away by moving them closer to the transmitter. • Check that the Repeater is set to the correct radio key. For further information see Function/Radio key.
The receivers in the system transmit signals for no reason.	<ul style="list-style-type: none"> • Change the Radio Key on all units in the system. For further information see Function/Radio key

For further information about the product in English, refer to the Appendix.

1. LED. Combined transmission and back-up battery indicator. For further information see Indicators and Signals
2. Test button/pushbutton
3. Exterior External Trigger.
4. Signal switch
5. Radio key switch
6. Connector for power supply unit
7. Battery cover
8. Wall mount bracket

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Appendix - further information

Connection

Connecting an Exterior External Trigger


The BE1510 Bellman Visit 868 Repeater has the following methods of activation:

- various accessories can be connected via the Exterior External Trigger socket (3), see below.
- with pushbutton

For more information about the indications provided by the receivers in the Bellman Visit 868 System, refer to Appendix/Further information/Settings below.

An Exterior External Trigger on the BE1510 Bellman Visit 868 Repeater can be connected to the vibrator output on several of Bellman & Symfon AB's products or to other equipment that produces a voltage or contact for activation. In addition, the following accessories can be used:

- magnetic switch BE9023 (accessory)
- door entering contact mat BE9024 (accessory)
- a standard pushbutton

1. Entry monitor
2. Magnetic switch
3. Pushbutton
4. Door phone
5. 3.5 mm jack plug
6. Exterior External Trigger

For further information see this appendix under the headline Technical information/Activation.

Settings

No adjustments are required for normal use. The relevant descriptions are provided below, if you wish to change a setting for some reason.

Radio key

In order to use several Bellman Visit 868 Systems close to one another without interference, different Radio Keys can be set on the different systems. All Bellman Visit 868 System units are supplied from the factory tuned to the same Radio Key, channel 0. This means that all Radio Key Switches on the transmitters are set to the OFF position.

- To alter the Radio Key, move the Radio Key Switches (5) to the desired positions.

Please note: all Bellman Visit 868 products within a system must be tuned to the same Radio Key in order to operate as a group.

Signal pattern

A Signal Pattern is the name for the way in which a receiver in the Bellman Visit 868 System indicates activation. Changing the transmitters' Signal Switch changes the Signal Pattern which the receivers display when the transmitter is activated. Follow the instructions under Appendix/Further information/Settings/Setting the activation pattern to make your choice.

The following signal patterns are available for the Bellman Visit 868 Repeater:

Type	LED pattern	Sound	Vibration	Flash
Green 1	Green is constantly lit	1 x ding dong, low-frequency tone	Separate	Yes
Green 2	Green blinks in sequences of two	2 x ding dong, low-frequency tone	Separate	Yes
Green 3	Green blinks in sequences of three	1 x ding dong, high-frequency tone	Separate	Yes
Green 4	Green blinks constantly	2 x ding dong, high-frequency tone	Separate	Yes
Yellow 1	Yellow is constantly lit	1 x ring, low-frequency tone	Short	Yes
Yellow 4	Yellow blinks constantly	2 x ring ring, high-frequency tone	Short	Yes
Orange 1	Orange is constantly lit	Baby	Rapid	Yes

More information about the signal patterns follow on the next side

Signal switch	Exterior external trigger	Pushbutton
	Green 1	Green 1
	Green 2	Green 2
	Green 3	Green 3
	Green 4	Green 4
	Orange 1	Green 4
	Orange 2	Green 3
	Orange 3	Green 2
	Orange 4	Green 1

Type	LED pattern	Sound	Vibration	Flash
Orange 2	Orange blinks in sequences of two	Baby	Rapid	Yes
Orange 3	Orange blinks in sequences of three	Baby	Rapid	Yes
Orange 4	Orange blinks constantly	Baby	Rapid	Yes
VMA	Red and Orange constantly blink alternately	VMA constant	Long	Yes
Fire alarm	Red blinks constantly	Fire alarm constant	Long	Yes

Setting the Activation Pattern

The Bellman Visit 868 Repeater can be activated in two different ways.

On delivery the Repeater sends the same signal regardless of the input causing the activation.

By changing the Signal Switches (4), different inputs on the same Repeater can transmit different Signal Patterns.

In this way the same Repeater can be used to forward signals from the Bellman Visit 868 System at the same time as it can transmit a Door Signal when activated via the Exterior External Trigger (3) and a Baby Cry signal when activated via the Test Button (2).

The following figure shows the standard settings that are available:

Testing

It is easy to test the Bellman Visit 868 Repeater. If the Repeater does not work as described below, you can check further under Troubleshooting/Troubleshooting guide.



How to test

A transmitter in the Bellman Visit 868 System which is set to the same radio key as the Repeater is required to test the radio reception on the BE1510 Bellman Visit 868 Repeater.

- Press the transmitter test button.
- The Bellman Visit 868 Repeater will indicate by blinking with the LED (1).

To test the Transmitter:

- Make sure that the Repeater is correctly connected.
- Activate the Repeater via the Test Button (2) or via the Exterior External Trigger (3).
- The Repeater should now light up LED (1) to show that it has been activated and that it is transmitting a signal to the receivers in the Bellman Visit 868 System.
- The receivers in the system will indicate an alarm according to the way that the Repeater has been set up using the Signal Switches (4).

Signal switch	Exterior external trigger	Pushbutton
	Green 1	Orange 1
	Green 2	Orange 2
	Green 3	Orange 3
	Green 4	Orange 4
	Green 1	Yellow 1
	Green 4	Yellow 4
	VMA	VMA
	Fire alarm	Fire alarm

Problem	Solution
The LED (1) blinks green when the Repeater is activated, but the receivers are not signalling.	<ul style="list-style-type: none"> • Check the battery in the receiver. • Check that all units in the Bellman Visit 868 System are set to the same radio key. For further information see Appendix/Further information/Settings/Radio key. • Check that the receivers are not placed too far away by moving them closer to the Repeater.
The receiver emits a strange alarm when activated by the Repeater.	<ul style="list-style-type: none"> • Set the Signal Switch to 0000.
The receivers in the system transmit signals for no reason.	<ul style="list-style-type: none"> • Change the Radio Key on all units in the system. For further information see Function/Radio key

Technical information

Power supply

Mains power: 6 V DC/800 mA with power supply unit BE9018 (Europe) or BE9083 (United Kingdom).

Back-up battery: 9 V 6LR61 alkaline
9 V 6F22 lithium

Power consumption: Active: 40 mA
Idle position: < 10 μ A

Troubleshooting

You can carry out a number of checks yourself before sending a product for repair.

Troubleshooting guide

Problem	Solution
Nothing happens when the transmitter is activated with the test button.	<ul style="list-style-type: none"> • Check that the power supply unit is connected correctly. • If the power supply unit is not connected, change the back-up battery. Only use an alkaline 6LR61 or a lithium 6F22 type battery. Make sure that the power supply unit is connected correctly. • Check that all the connections are correct.
The LED (1) blinks yellow when the Repeater is activated.	<ul style="list-style-type: none"> • Change the back-up battery. Only use an alkaline 6LR61 or a lithium 6F22 type battery. Make sure that the power supply unit is connected correctly.

The troubleshooting guide continues on the next page.

Activation

Via test button

- Exterior External Trigger:** 3.5 mm stereo (mono provides a connection) jack plug (3)
- Connection:** Between the inner and outer pins of the mono type 3.5 mm jack plug or between the middle/inner and outer pins of the stereo type 3.5 mm jack plug, see diagram.
- DC:** 2 to 30 V between the inner pin and middle pin on the stereo type 3.5 mm jack plug, see diagram.
- AC:** 3 to 24 V RMS 5 -150 Hz between the inner pin and middle pin on the stereo type 3,5 mm jack plug, see diagram.

Additional information

For indoor use only

Size B x H x D: 60 x 98 x 30 mm**Weight:** With back-up battery: 110 g
Without back-up battery: 70 g**Colour:** White with grey pushbutton**Accessories**BE9024 Door Entering Contact Mat
BE9023 Magnetic Switch**Radio function****Radio frequency:** 868.3 MHz

Number of Radio Keys: 64 Radio Keys as standard. Special software can be used to increase these to 256 Radio Keys in increments of 64 per software purchase. Contact the nearest supplier for further information.

Coverage: The normal coverage between a transmitter and receiver in the BellmanVisit 868 System is approximately 200 metres with a clear line of sight. Coverage is reduced if walls and large objects screen off the signals. Any thick walls constructed of reinforced concrete will greatly affect coverage. The system may also be affected by radio transmitters such as TV transmitters, computers, mobile phones, etc. This means that a unit may work perfectly in one part of the room but not at all in another.