Maxi

Personal amplifier

DE

For use in all EU countries

Hereby, Bellman & Symfon Group AB, declares that this BE2020 is in compliance with the essential requirements and other relevant provisions of: EMC directives: 2014/30/EU. MDD directives: 93/42/EEC. RoHS directive: 2015/863/EU. EN 60601-1-2, Class B Passed.

Correct use and disposal of batteries:

Replace only with the same or equivalent type recommended by the manufacturer. Please dispose of old, defective batteries in an environmentally friendly manner in accordance with the relevant legislation.

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Appendix

Connection

Normally when you use Bellman Audio Maxi digital communication aid, you place it on a table, hold it in your hand, hang it round your neck using a neck strap (accessory) or attach it to an article of clothing with the clip supplied (18).

You can connect BE9122 Bellman Audio Headphones, BE9159/BE9161 Bellman Audio Neck Loop, BE9403 Bellman Audio Stetoclips, BE9124 Bellman Audio Ear Phones or another product recommended by Bellman & Symfon to the Headphone Output (12).

You can connect a BE9126 Bellman Audio Cable Kit to the External Sound Source Input (1), which you can connect to the required sound source or a BE9127 Bellman Audio External Microphone. It has a 5 metre long cable so that it can, for example, be positioned next to the TV or someone who is speaking.

Press on the On/Off button  (5) to start Maxi. While Maxi is starting up, the Microphone LED (3) and the Telecoil LED (4) blinks. When the Microphone LED (3) comes on, Maxi is ready to use and you will hear sound from the internal microphone in the headphones.

Settings

Bellman Audio Maxi has more setting options than those found on the front. You can adjust and adapt Maxi to your individual requirements using the Setting button (15).

Balance

The balance between the right and left channels can be set, e.g. for headphones. This setting can be used to compensate for a hearing loss in one ear.

The standard setting is centred balance.

All settings are carried out in sequence, see below under Adjusting settings.

Basic attenuation

The basic attenuation can be set, for example, for different headphones. What you actually set is the attenuation you require at the maximum power output you want to have (which is therefore 0 dB attenuation). This setting is useful when you want to adjust the maximum volume for a pair of headphones or in-ear phones with a higher sensitivity than the BE9122 Bellman Audio Headphones.

The standard setting is 0 dB attenuation.


All settings are carried out in sequence, see below under Adjusting settings.

Blocking the T position

The Telecoil Selection Button can be blocked if you know that the T position will not be used. This setting is for those who will never use the Telecoil and want to avoid selecting this position by mistake.

The standard setting is for blocking of the T position not to be activated.



All settings are carried out in sequence, see below under Adjusting settings.







Adjusting settings

It is easy to adjust settings.

- Open the battery cover (17) by pressing the battery cover lock (16) in the direction of the arrow.
- Press the Setting button (15) for about 5 seconds. It is easiest to do this with a paper clip, see figure. One of the LEDs in the Volume and Tone LED Indicator (8) now starts to blink
- Adjust the balance between the right and left channels using the volume







control  (9).  moves the sound towards the right channel and  moves the sound towards the left channel. The current setting is indicated by a flashing dot on the Volume and Tone LED Indicator (8) which moves towards the right or left.

When the desired value has been set, give a short press on the Setting button (15) and Maxi will move on to setting the maximum volume.

- Adjust the maximum volume using the volume control  (9).  increases the maximum volume and  reduces the maximum volume. The current attenuation setting relative to the maximum volume is displayed via the Volume and Tone LED Indicator (8). No attenuation is shown by only one LED being visible at the right. The greater the attenuation the more LEDs come on from right to left. Each increment and hence each lit LED corresponds to a reduction of 3 dB. When the desired value has been set, give a short press on the Setting button (15) and Maxi will move on to the setting for locking the T position.
- Now set the use of the telecoil to locked or unlocked using the volume control  (9).  enables use of the telecoil and  blocks use of the telecoil. When the telecoil is activated, the Telecoil LED (4) blinks and when the function is blocked the Telecoil LED (4) is off.
- When the desired value has been set, there are three options:
 - o To start again: give a short press on the Setting button (15) and Maxi starts again with setting the balance above.
 - o To save: press the Setting button (15) for about 5 seconds and Maxi saves the settings that have been adjusted and returns to normal mode.
 - o To change your mind: press the On/Off button  (5) and Maxi will disregard all the settings that have been adjusted and return to normal mode.

Tone control lock

To prevent that the tone control is accidentally changed by the user it is possible to lock the tone control.

1. While the Maxi is off press and hold the  button (10). Switch on the Maxi, by pressing the On/Off button  (5). Keep the  button held down until the Maxi's lights (3, 4) have stopped flashing.
2. One of the LED in the Volume and Tone LED Indicator (8) will be flashing, this will be the set tone level. Adjust the preset tone level by using the Tone control   (10). If all the LED's on the Volume and Tone LED Indicator (8), except the blinking one, are lit at the same time, it means that the Tone control is activated while if the rest of the LED's are off the Tone control is deactivated.
3. Press the  button to activate the Tone control or press the  button to deactivate the Tone control.
4. Switch the Maxi off using the On/Off button  (5) to save these settings.

Testing

It is easy to test the Bellman Audio Maxi digital communication aid for yourself. If Maxi is not working as described below, you can carry out further troubleshooting as instructed in the section Troubleshooting/Troubleshooting guide.

How to test

The internal microphone and telecoil can be tested by connecting BE9122 Bellman Audio Headphones to the Headphone output (12) and then selecting the Internal microphone (2) using the Microphone Selection Button  (6) or the Telecoil Selection Button  (7).

The External sound source input (1) can be tested by connecting a CD player or a radio to the input via an audio cable. Switching to the External signal source input (1) takes place automatically when the connector is connected to the External sound source input (1) on Maxi.

Troubleshooting

You can carry out a number of checks yourself before sending a product for repair.

Troubleshooting guide

Symptom	Solution
Nothing happens when you try to start Maxi using the On/Off button.	<ul style="list-style-type: none"> • Check that the batteries have been inserted the right way round. • Change the batteries. Only use AA alkaline type batteries.
No sound can be heard in the headphones when Maxi is switched on.	<ul style="list-style-type: none"> • Check that the correct signal source has been selected, e.g. internal microphone (M position). • Check that the headphone plug is securely inserted in the headphone output. • Connect the headphones to another sound source (with a 3.5 mm headphone output) to test the headphones.
A high-pitched noise is heard in the headphones.	<ul style="list-style-type: none"> • Lower the volume or increase the distance between Maxi and the headphones. • Direct the microphone away from you.

Technical information

Power supply

Battery power:	Two 1.5 V AA alkaline type batteries
Operating time:	100 - 150 hours depending on the load and sound pressure
Power consumption:	17 mA (120 mA max)

Input signals

Microphone:	Built-in microphone (omnidirectional)
Telecoil:	Built-in telecoil (angled at 45°)
Ext. sound source input:	2.5 mm stereo jack plug
Ext. in sensitivity (max input lvl.):	+6 dBV

Output signals

Headphone socket:	3.5 mm stereo jack plug
Output level with BE9122:	117.67 dB @ 1 kHz (SPL90) 119.52 dB (HF Ave. SPL90), 37.01 dB
Full-on gain:	29.21 dB (SPL70)
Ref. test gain:	148.7 dB @ 1 kHz (SPL90)
Output level with BE9403:	52.6 dB (HF Ave. SPL90)
Ref. test gain:	71 dB (SPL70)
Distortion:	0.559% THD (electrical)
SNR:	82 dB
Frequency range:	40 Hz – 10 kHz
Output power:	Max 125 mW @ 16 ohms

Features

Volume control:	9 steps (5 dB/step)
Maximum volume settings:	0-21 dB programmable attenuation (in settings mode only and in steps of 3 dB)
R/L-channel balance setting:	0-20 dB programmable attenuation (in steps of 5 dB)
Tone control:	+/- 10 dB (5 pre-set steps, push-buttons) <i>Tone, high:</i> +10dB@3.15 kHz <i>Tone, mid:</i> -6dB@250Hz <i>Tone, low:</i> -10dB@3.15 kHz No bass boost
Power switch:	on-off toggle button
Mic and T-mode:	Separate mode select buttons
Ext. mic / Aux in accessories:	Automatically selects Aux in or ext. mic. when plugged in
MicroSet™ accessory:	Automatically selects MicroSet™ accessory when plugged in
Dynamic compressor:	10 channels Dynamic gain range: 0-35 dB Compression: 3:1 Expansion: 1:1.5 Attack time: 25 ms Release time: 250 ms

Dynamic noise reduction:
Feedback cancellation:
Digital signal processing:

10 channels Adaptive frequency domain noise reduction
Adaptive time domain feedback cancellation
19.948 kHz sampling frequency 16-bit resolution in stereo

Additional information

For indoor use only.	
Dimensions WxHxD:	45 x 140 x 27 mm
Weight (without clip):	With battery: 133 g Without battery: 83 g
Colour:	Grey with white front panel and grey buttons.
Operating temperature:	0°– 35° C, 32°– 95° F.
Transport and storage temperature:	-10°– 50° C.
Relative humidity:	15% – 90% (non condensing).

Regulatory requirements

FCC SDoC: FCC Part 15 Subpart B,
CE, RoHS, WEEE, RCM

Accessories

BE9122	Bellman Audio Headphones
BE9159/BE9161	Bellman Audio Neck Loop
BE9403	Bellman Audio Stetoclips
BE9124	Bellman Audio Ear Phones
BE9126	Bellman Audio Cable Kit
BE9127	Bellman Audio External Microphone

Measurements

Frequency response (acoustical):