
# Visit Flash Receiver

Model BE1441/1442/1443

CZ

DE

DK

ES

FI

FR

GB

GR

HU

IT

NL

NO

PL

PT

SE

SRB

APP

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GB

### Settings

**Please note:** No adjustments are required for normal use. To change a settings, follow the steps below.

#### Changing the radio key

If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference, you need to change the radio key on all units. The radio key switches are located on the transmitters.

#### Here is how you change the radio key:

- 1 Open the transmitter cover and move any radio key switch to the up = on position to change the radio key. See **Changing the radio key** for the relevant transmitter.
- 2 Press and hold the mute/test button on the top of the receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- 3 Press the test button/s on the transmitter within 30 seconds to send the new radio key.
- 4 All Visit LEDs on the receiver blink 5 times to show that the radio key has been changed. It then returns to normal mode.

**Please note:** All Visit units must be set to the same radio key in order to operate as a group.

#### Changing the signal pattern

When a transmitter is activated, the flash receiver lights up a LED, flashes and the bed shaker starts to vibrate with a certain pace. This is called signal pattern. The transmitters determine the pattern. Changing the transmitters' signal switch changes the signal pattern on the flash receiver.

The following signal patterns are available for the Visit system:  
(Please note that the flash receiver does not emit any sound.)

#### 868/315MHz Europe, Asia


| Type       | LED pattern                                 | Sound                              | Vibration | Flash |
|------------|---|------------------------------------|-----------|-------|
| Green 1    | Green is constantly lit                     | 1 x ding dong, low-frequency tone  | Separate  | Yes   |
| Green 2    | Green blinks in sequences of two            | 2 x ding dong, low-frequency tone  | Separate  | Yes   |
| Green 3    | Green blinks in sequences of three          | 1 x ding dong, high-frequency tone | Separate  | Yes   |
| Green 4    | Green blinks constantly                     | 2 x ding dong, high-frequency tone | Separate  | Yes   |
| Yellow 1   | Yellow is constantly lit                    | 1 x ring, low-frequency tone       | Short     | Yes   |
| Yellow 2   | Yellow blinks in sequences of two           | 2 x ring ring, low-frequency tone  | Short     | Yes   |
| Yellow 3   | Yellow blinks in sequences of three         | 1 x ring, low-frequency tone       | Short     | Yes   |
| Yellow 4   | Yellow blinks constantly                    | 2 x ring ring, high-frequency tone | Short     | Yes   |
| Orange 1   | Orange is constantly lit                    | Baby                               | Rapid     | Yes   |
| Orange 2   | Orange blinks in sequences of two           | Baby                               | Rapid     | Yes   |
| Orange 3   | Orange blinks in sequences of three         | Baby                               | Rapid     | Yes   |
| Orange 4   | Orange blinks constantly                    | Baby                               | Rapid     | Yes   |
| VMA        | Red and Orange constantly blink alternately | VMA constant                       | Long      | Yes   |
| Fire alarm | Red blinks constantly                       | Fire alarm constant                | Long      | Yes   |

## 433MHz North America, Australia

| Type       | LED pattern                                 | Sound                      | Vibration | Flash |
|------------|---|----------------------------|-----------|-------|
| Green 1    | Green is constantly lit                     | 1 x Ding dong, Bass tone   | Long      | Yes   |
| Green 2    | Green is constantly lit                     | 1 x Ring, Treble tone      | Short     | Yes   |
| Green 3    | Green is constantly lit                     | 1 x Ding dong, Bass tone   | Long      | Yes   |
| Yellow 1   | Yellow is constantly lit                    | 1 x Ding dong, Treble tone | Long      | Yes   |
| Yellow 2   | Yellow is constantly lit                    | 1 x Ring, Bass tone        | Short     | Yes   |
| Yellow 3   | Yellow is constantly lit                    | 1 x Ding dong, Treble tone | Long      | Yes   |
| Orange 1   | Orange is constantly lit                    | 2 x Ding dong, Bass tone   | Long      | Yes   |
| Orange 2   | Orange is constantly lit                    | 2 x Ring, Bass tone        | Short     | Yes   |
| Red 1      | Red is constantly lit                       | 2 x Ding dong, Treble tone | Long      | Yes   |
| Red 2      | Red is constantly lit                       | 2 x Ring, Treble tone      | Short     | Yes   |
| Red 3      | Red is constantly lit                       | 2 x Ding dong, Treble tone | Long      | Yes   |
| VMA        | Red and Orange constantly blink alternately | VMA                        | Constant  | Yes   |
| Fire alarm | Red blinks constantly                       | Fire                       | Constant  | Yes   |

## Battery backup

BE1442



**CAUTION**

RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE.

DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS.


**Getting started**

- Remove the battery slip on the bottom of the unit before use.
- Please charge the unit for 24 hours to ensure full backup time.

**Replacing batteries  
(Done by the operator)**


- Use only 4x1.2VNiMH 600AAA/900AAA Lexel rechargeable batteries to ensure full backup time.

## Battery backup

### Getting started

- Remove the battery slip on the bottom of the unit before use.
- Please charge the unit for 24 hours to ensure full backup time.


### Replacing batteries (Done by the operator)

- Remove the screw on the bottom of the unit to exchange the batteries.

#### CAUTION

RISK OF EXPLOSION IF BATTERY IS REPLACED BY AN INCORRECT TYPE.  
DISPOSE OF USED BATTERIES ACCORDING TO THE INSTRUCTIONS.


- Use only 4x1.2VNiMH 900AAA Lexel rechargeable batteries to ensure full backup time.

## BE1443


## Smoke detector supervision

Pairing a smoke detector with the Flash receiver

**1**

**2**


- Press (9) and hold the top button for 5 seconds.
- The flash receiver indicates that it is ready by switching between the green and the yellow LED light.
- Release the top button.

- Press and hold the test button on the smoke detector for more than 3 seconds.
- Release the button.
- The flash receiver will indicate that pairing was successful by switching between the orange, green, yellow and red LED lights three times ☺💡🔥.

## BE1443

**i**

- This menu shows the position of the smoke detector.
- Attach the corresponding label (provided with each smoke detector) on the smoke detector unit so you can easily identify it later.

## Flash receiver error identification

**BE1443**

If a paired smoke detector is faulty, the red LED  light on the Flash receiver will start to blink. By looking at the power LED  and the smoke detector LED  you can easily identify what is wrong. If you are using multiple smoke detectors, just enter the smoke detector menu to find the faulty unit.

| Error            | Indicator  |
|------------------|--|
|                  | Power LED  Smoke detector LED  |
| Low battery      | • RED, blinking<br>• AMBER   |
| Lost connection  | • GREEN<br>• AMBER, blinking   |
| Hardware failure | • OFF<br>• AMBER, blinking   |


## Trouble shooting

Most problems can be solved quickly by following the advice below.


| If                                       | Try this   |
|--|--|
| The smoke detector indicates low battery | • Replace the battery using the correct batteries stated at the smoke alarm label.   |
| The connection is lost                   | • Check if the battery in the smoke detector is still OK.<br>• Move the units to their original position and verify that the error message disappears. |
| There is a hardware failure              | • Clean the smoke detector according to the smoke detector manual.<br>• If it is still faulty, replace the smoke detector.                             |

**BE1443**

## Removing one smoke detector from the system


- Press the Smoke detector Menu button  on the backside of the unit.
- The Smoke detector LED  will blink green to indicate you are entering the menu system.
- Press the Smoke detector Menu button  as many times as needed to reach the smoke detector you would like to remove from the system.


- Press and hold the Smoke detector Menu button for more than 5 seconds.
- Release the button.
- The flash receiver will indicate the removal of the smoke detector by switching between the orange, green, yellow and red LED lights three times.

## Removing all smoke detectors from the system

- Press the Smoke detector Menu button  on the backside of the unit.
- The Smoke detector LED  will blink green to indicate you are entering the menu system.

- To remove all smoke detectors in the system, press and hold the Smoke detector Menu button  for more than 5 seconds.
- Release the button.
- The flash receiver will indicate the removal of the smoke detector by switching between the orange, green, yellow and red LED lights three times.

## Advanced programming

By using advanced programming, you can customize the signal pattern from a specific transmitter and event, displaying the LED colour and vibration pattern of your choice. The advanced programming overrides the radio key and pairs the units via the serial number. Please note that smoke alarms cannot be programmed for safety reasons.

**Note:** The transmitter must be activated as it is intended to be used in the system to generate the right signal.

This means that you can't always use the transmitter test button (see Default signal pattern for the relevant transmitter).

### Here is how you program the receiver:

- 1 Press and hold the mute/test button on the receiver. The green and yellow Visit LEDs will start to blink alternately. While still holding down the button, activate the desired transmitter as intended. Release the button.
- 2 Scroll through the different **Visit LED options** by pressing the mute/test button on the receiver. Select the desired Visit LED colour by holding down the mute/test button until the power LED goes out and lights up again.
- 3 Scroll through the different **vibration options** by pressing the test button on the receiver (bed shaker required). Select the desired vibration pattern by holding down the mute/test button until the power LED goes out and lights up again.
- 4 The receiver will now show the new Visit LED colour and vibration pattern. Press the mute/test button briefly to end the demonstration. After a short while, it will return to normal mode.

### Deleting the advanced programming

Follow the procedure below to delete the advanced programming.

- 1 Hold down the mute/test button on the receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- 2 Press the mute/test button on the receiver 3 times in quick succession.
- 3 All Visit LEDs will light up for ~2 seconds to show that it has been deleted.

## Testing

It is easy to test the flash receiver. If it does not work as described below, see the section **Troubleshooting**.

### Testing the system

A Visit transmitter set to the same radio key as the flash receiver is required to test the system. You will need a bed shaker accessory to test the vibration.

- 1 Press the transmitter test button.
- 2 The flash receiver will:
  - Start to flash
  - Light up the corresponding Visit LED
  - If a bed shaker accessory is connected, it will start to vibrate

### Testing the landline telephone connection

- 1 Connect the flash receiver to an analogue telephone socket via the telephone Input. You will need a telephone flex cable accessory and an adapter plug accessory.
  - 2 Use for instance your mobile phone to call your landline telephone. The flash receiver will start to flashing and light up the yellow Visit LED. If a bed shaker accessory is connected, it will start to vibrate.
- To repeat the last indication, press the flash receiver test button.

## Troubleshooting

Most problems with the flash receiver can be solved quickly by following the advice below.

| If   | Try this  |
|--|---|
| Nothing happens.                                       | <ul style="list-style-type: none"> <li>▪ Check that the power supply unit is connected correctly. The power LED should light up green.</li> <li>▪ Check that there is current in the wall socket.</li> </ul>  |
| The receiver is not activated.                         | <ul style="list-style-type: none"> <li>▪ Check the battery in the transmitters.</li> <li>▪ Check that the receiver is not placed too far away from the transmitters by moving it closer to the transmitters.</li> <li>▪ Check that the receiver is set to the correct radio key. For further information see <b>Radio key</b>.</li> </ul> |
| The receiver signals when no transmitter is activated. | <ul style="list-style-type: none"> <li>▪ Change the radio key on all units in the system. For further information see <b>Radio key</b>.</li> </ul>  |

## Technical information

### Power supply

|             |  |
|-------------|--|
| Mains power | 7.5 V DC / 1500 mA<br><b>Note:</b> Only use the supplied power adapter.<br>This product is also designed for IT power distribution system with phase to phase voltage of up to 230V. The socket-outlet shall be installed near the equipment and shall be easily accessible. |
|-------------|--|

|                |  |
|----------------|--|
| Battery backup | For BE1441: not applicable<br>For BE1442: 4x1.2V Lexel 600/900mAh<br>For BE1443: 4x1.2V Lexel 900mAh |
|----------------|--|

Please dispose used batteries at municipal collecting points or dispense them to local stores free of charge.

### Radio function

|                   |   |
|-------------------|---|
| Power consumption | Active current: 900mA idle current: 10mA    |
| Radio frequency   | 314.91/433.92/868.3 MHz depending on region |

|                      |   |
|----------------------|---|
| Number of Radio Keys | 64 Radio Keys as standard. Special software can be used to increase these to 256 Radio Keys in increments of 64 per software purchase. Contact the nearest supplier for further information.  |
| Coverage             | The normal coverage between a transmitter and receiver in the Bellman Visit System is up to 250 metres (868MHz) / 80m (433MHz) / 50m (314.91MHz) with clear line of sight. Coverage is reduced if walls and large objects screen off the signals. Any thick walls constructed of reinforced concrete will greatly affect coverage. The system may also be affected by radio transmitters such as TV transmitters, computers, mobile phones, etc. This means that a unit may work perfectly in one part of the room but not at all in another. |

### Additional information

|                      |  |
|----------------------|--|
| Programming pin (16) | For service personal only (Output: 3.3V/max. 150mA)<br><b>Caution:</b> Do not use, unit will be damaged when used. |
|----------------------|--|

### Activation

|                            |   |
|----------------------------|---|
| Radio                      | Visit system  |
| Via analogue telephone nw. | 26 - 120 V RMS, 15 - 100 Hz.  |
| BE1441/BE1442              | Via the Mobile phone detector, see the accessories section (Output: 3.3V/max 150mA) |

**Warning!** The telephone port of the Flash receiver must be connected using 26 AWG or greater telecommunications wire.

### Output signals

|                |  |
|----------------|--|
| Output signals | Built-in flash light signal ~30 Candela<br><b>Warning!</b> Flashes can trigger epileptic seizures. |
|----------------|--|

|                            |  |
|----------------------------|--|
| Vibrator power             | 2.0 – 4.0 VDC (max. 500mA)   |
| Environmental requirements | Operating temperature: 0°C – 35°C<br>Transport and storage temperature: -10°C – 50°C<br>Relative humidity: 15% - 90%, non condensing |
| Dimensions Ø x H           | BE1441: 70 x 140 mm<br>BE1442/BE1443: 78 x 155 mm  |

|        |                               |
|--------|-------------------------------|
| Weight | 308 g, 10.8 oz (three models) |
| Colour | White                         |

|                      |   |
|----------------------|---|
| Optional accessories | BE9075 Wall mount bracket<br>BE1270 Bed shaker<br>BE9105 Telephone flex cable<br>BE9250 Mobile phone sensor<br>Adapter plug for the appropriate country |
|----------------------|---|

## FCC compliance statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**Warning:** Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**Note:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

## Statement for Industry Canada

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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