Visit Polsontvanger & Lader

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Settings

No adjustments are required for normal use. The relevant descriptions are provided below, if you wish to change a setting for some reason.

Radio key

In order to use several Visit systems close to one another without interference, different radio keys can be set on the different systems. All Visit system units are supplied from the factory tuned to the same radio key, channel **0**. This means that all radio key switches on the transmitters are set to the **OFF** position.

Changing radio key

A Visit transmitter is also required to alter the radio key on the wrist receiver. Proceed as follows to change the radio key:

- 1 Set the transmitter to the desired radio key by altering its radio key switch to the desired radio key. See the transmitter user manual for further information.
- 2 Hold down the function button (1) on the wrist receiver until the green  (3) and yellow  (4) LEDs blink alternately. Release the button and the wrist receiver will now be in programming mode for about 30 seconds.
- 3 Press the test button on the transmitter. The LEDs     (2, 3, 4 and 5) will now blink five times in quick succession to indicate a successful change of radio key.
- 5 After changing the radio key, the wrist receiver will automatically return to normal mode again.

Please note: all Visit products within a system must be tuned to the same radio key in order to operate as a group.

Type	LED pattern	Vibration pattern
Green 1	Green constant light	Slow
Green 2	2 x green blinks	Slow
Green 3	3 x green blinks	Slow
Green 4	Blinking in green	Slow
Yellow 1	Yellow constant light	Medium
Yellow 2	2 x yellow blinks	Medium
Yellow 3	3 x yellow blinks	Medium
Yellow 4	Blinking in yellow	Medium
Orange 1	Orange constant light	Fast
Orange 2	2 x orange blinks	Fast
Orange 3	3 x orange blinks	Fast
Orange 4	Blinking in orange	Fast
VMA	Blinking in red and orange	Long
Fire alarm	Blinking in red	Long
Fire alarm low bat.	Blinking in red every 5s	None
Fire alarm flat bat.	Blinking in red every 5s	One short

● LED blink ○ LED off ●→ Constant LED light ■ Vibration, Pause

Advanced programming

Advanced programming provides additional options for those who wish to make special modifications to the wrist receiver.

The idea is that it should be possible to select a completely unique signal pattern which is linked to activation from a specific input on a specific transmitter. The function works regardless of the radio key settings on the units that are programmed. Please note that, for safety reasons, the function does not work with the Visit smoke alarm transmitter.

By using advanced programming of the wrist receiver, it can be adjusted so that its signal pattern corresponds exactly to what is required. In other words, an entirely individual signal pattern can be programmed, such as displaying a permanently lit orange LED and a constant vibration.

In order to adjust the setting, the transmitter to which the wrist receiver should be adapted must be available. The transmitter must also be connected so that it can be activated in the way in which it is intended to be used.

Proceed as follows

- 1 Hold down the function button (1) on the wrist receiver until the green  (3) and yellow  (4) LEDs blink alternately. The wrist receiver will now be in programming mode for about 30 seconds.
- 2 Hold down the function button (1) at the same time as the relevant transmitter is activated in the way in which it is intended to be used. The yellow LED  (6) will light up. Note that all inputs are individual, so that it is not possible to use the test button on the telephone transmitter whose telephone input will indicate the relevant pattern.
- 3 Scroll through the different LED options by a short press on the function button (1). Select the relevant indication by holding down the function button (1) until the LED  (6) goes out and starts to shine with a constant yellow light again.
- 4 Scroll through the different vibration options by a short press on the function button (1). Select the relevant indication by holding down the function button (1) until the LED  (6) goes out and starts to shine with a constant yellow light again.
- 5 The wrist receiver will now show the indication method programmed. End the display with a short press of the function button (1).
- 6 After a short while the wrist receiver will automatically return to normal mode.

Resetting advanced programming

It is quite easy to reset the wrist receiver if it needs to be reset after it has been programmed

using advanced programming.

- 1 Hold down the function button (1) on the wrist receiver until the green  (3) and yellow  (4) LEDs blink alternately. The wrist receiver will now be in programming mode for about 30 seconds.
- 2 Press the function button (1) three times in quick succession.
- 3 All LEDs     (2-5) remain constantly on for a few seconds.
- 4 All the advanced programming has now been deleted and the wrist receiver will automatically return to normal mode.

Testing the wrist receiver

It is easy to test the wrist receiver. If it does not work as described below, you can check further below in the **Troubleshooting** guide in the appendix. A transmitter which is set to the same radio key as the wrist receiver is required in order to test the radio reception.

Press the transmitter test button. The wrist receiver will give the following indications:

- The LEDs     (2-5) which the transmitter has been set to indicate with ,or which have been programmed with advanced programming, will be switched on.
- It will vibrate as the transmitter is set to indicate, or in the way the wrist receiver has been programmed with advanced programming.

Troubleshooting

If	Try this
<p>Nothing happens</p>	<ul style="list-style-type: none"> ▪ Check that the power supply unit is connected correctly and that the green LED (10) on the charger lights up. ▪ Charge or renew the battery in the wrist receiver. Only use a rechargeable VARTA V40H NiMH battery. If the battery is completely flat, you can still use the wrist receiver by placing it in the charger. The wrist receiver needs to be charged for about two hours before it can be used separately. ▪ Check that the wrist receiver has been correctly placed in the charger and that the LED (12) lights up in green.
<p>It is not possible to charge the wrist receiver</p>	<ul style="list-style-type: none"> ▪ Check that the power supply unit is connected correctly and that the green LEDs (10, 12) on the charger light up. ▪ Charge or change the batteries in the charger. Only use rechargeable TMK AAA60CT NiMH batteries.
<p>Nothing happens when the power supply unit is disconnected</p>	<ul style="list-style-type: none"> ▪ Charge or change the batteries. Only use rechargeable TMK AAA60CT NiMH batteries.
<p>The LED  (6) on the wrist receiver blinks yellow</p>	<ul style="list-style-type: none"> ▪ Charge or renew the battery in the wrist receiver. Only use a rechargeable VARTA V40H NiMH battery. ▪ If the battery is completely flat, you can use the wrist receiver by placing it in the charger. It needs to be charged for about 2 hours before it can be used separately.

If	Try this
<p>The bed shaker does not vibrate when the wrist receiver is activated</p>	<ul style="list-style-type: none"> ▪ Check that the bed shaker is correctly connected to the charger. ▪ Check that the wrist receiver has been correctly inserted in the charger and that the LED (12) lights up green.
<p>The wrist receiver does not pick up signals from the transmitters</p>	<ul style="list-style-type: none"> ▪ Check the batteries in the transmitters. ▪ Check that the wrist receiver is not placed too far away by moving it closer to the transmitter. ▪ Check that the wrist receiver is set to the same radio key as the other units in the relevant Visit system. For further information see Changing radio key in the appendix.
<p>The wrist receiver is emitting the wrong signals, regardless of how the transmitter is set</p>	<ul style="list-style-type: none"> ▪ Delete the advanced programming, see Resetting advanced programming in the appendix.
<p>The receiver signals when no transmitter is activated</p>	<ul style="list-style-type: none"> ▪ Change the radio key on all units in the relevant Visit system. There is probably another Visit system nearby using the same radio key.

Technical information

	BE1560 Visit wrist receiver	BE1570 Charger
Mains power	Via the charger, BE1570	7.5 V DC/1500 mA Only use the supplied power adapter
Battery power	1 x 1.2 V VARTA V40H rechargeable NiMH battery. The battery cover can be opened with a coin.	4 x 1.2 V/600 mA AAA NiMH TMK AAA60CT rechargeable batteries
	<p>Please note!</p>  Only use rechargeable VARTA V40H NiMH batteries!	<p>Please note!</p>  Only use rechargeable VARTA TMK AAA60CT NiMH batteries!
Operating time	Approx. 30 hours	-
Charging time	Approximately 8h. If the battery is completely flat, you can still use the wrist receiver by placing it in the charger. It needs to be charged for about 2h before it can be used separately.	Discharged backup batteries will take at least 24h to charge. The batteries are charged via the power supply unit. The charger should normally be permanently connected to mains power.
Power consumption	Active: 100 mA Idle position: 3 mA	Active: 650 mA Idle position: 70 mA
Radio frequency	868.3 MHz	No radio
No. of radio keys	64 radio keys as standard Up to 256 with special software	-
Coverage	~100m, clear line of sight	-
Radio activation	The Visit 868 system	-
Vibrator power	-	2.0 - 4.0 VDC
Flex length	-	1.7 m
Dimensions	38 x 12 x 49 mm (w x h x d)	95 x 100 x 117 mm (w x h x d)
Weight	27g Incl. battery	185g incl. battery
Accessories	BE1570 Charger BE1270 Bed shaker BE9140 Wrist belt tool	BE1270 Bed shaker BE9086 External trigger cable

Additional information

Operation, storage and transport conditions Operating temperature: 0°–35° C, 32°–95° F. Transport and storage temperature: -10°–50° C. Relative humidity: 15% – 90% (non condensing).

Hereby, Bellman & Symfon, declares that this wrist receiver and charger is in compliance with the essential requirements and other relevant provisions of directives: RED directive 2014/53/EU.

For use in all EU countries

Correct use and disposal of batteries
 Replace only with the same or equivalent type recommended by the manufacturer.
 Please dispose of old, defective batteries in an environmentally friendly manner in accordance with the relevant legislation.

Contact us

For a complete Declaration of Conformity please contact info@bellman.com or call +46 31 68 28 20

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Bellman & Symfon Europe AB
Bellman & Symfon Asia Ltd.