Important safety instructions

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Warning!

- This device has a very powerful amplifier and the sound level can cause discomfort or, in serious cases, endanger hearing if the product is used carelessly. Always turn the volume down to a low level before putting on the headphones!
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Maxi Pro BE2021, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Failure to follow these safety instructions could result in fire, electric shock, or other injury or damage to the product or other property.
- Do not use or store this product near any heat sources such as naked flames, radiators, ovens or other devices that produce heat.
- Protect the cables from any potential source of damage.

- Use only power adapters and accessories specified in this manual.
- Only refer a damaged or faulty product to qualified service personnel.
- Do not dismantle the product; there is a risk of electric shock. Tampering with or dismantling the product will void warranty.
- This product is designed for indoor use only.
- Do not expose the product to moisture.
- Protect the product from shocks during storage and transport.

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Battery charging and replacement

Warning! The Maxi Pro is powered by a rechargeable lithium-polymer battery. Do not expose the battery to fire or to direct sunlight. The battery should only be changed by Bellman & Symfon or an authorized service center. There is a risk of explosion if the battery is not replaced correctly. The battery should be charged fully before using the product for the first time. Only charge the battery in a temperature between 0°– 35°C, 32°– 95° F.

Using the product

Operating environment: 0°– 35° C with an atmospheric pressure of 700hpa - 1060hpa.
Transport and storage environment: -10°– 50° C with an atmospheric pressure of 700hpa - 1060hpa. Do not use or store the products near a heat source. If the product gets wet or is exposed to moisture, it should be replaced.

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Technical specifications

Intended Use

Personal sound amplification product that enhances speech to facilitate listening during conversations, telephone calls and television viewing.

Dimensions and Weight

Size: 140 × 27 × 45 mm, 5.5 × 1.1 × 1.8"
Weight: 127g, 4.3 oz.

Frequency and Coverage

Bluetooth frequency: 2402 – 2480 MHz
Coverage: Up to 25 m, 27.5 yd, free field.
Connectivity: Qualcomm® aptX™ LL

Power Supply

Input: 100 – 240V~, Output: DC 5V, 1.0A
Battery power: 2 × 3.7V 750mAh Li-on
Operating time: Up to 70 h
Power consumption: 18 mA

Connections

Headphones input: 3.5 mm mini-stereo jack
Charger input: 2.5 mm jack

Microphone

Built-in omnidirectional microphone

Output level with BE9229

135dB @ 1kHz (SPL90)
136dB (HF Ave. SPL90)
Full-on gain: 49dB @ 1kHz,
Ref. test gain: 49dB @ 1kHz (SPL70)

Configuration options

These devices can be configured with the following original accessories:

- BE9233 Headset with microphone, black
- BE9229 Earphones with mic. black
- BE9159 Neck loop, 90 cm cable, black
- BE9161 Neck loop, 15 cm cable, black
- BE9254 Single charger cord, white
- Power adapter, white
- USA plug, white
- AU plug, white

Output level with BE9233

122dB @ 1kHz (SPL90)
127dB (HF Ave. SPL90)
Full-on gain: 36dB @ 1kHz
Ref. test gain: 36dB @ 1kHz (SPL70)

Mic SNR

53dB @ 1kHz 1% THD (electrical)

Distortion

0.2% THD @ 1kHz – 6dB @ 1% THD (Mic electrical)

Bluetooth SNR

85dB @ 1kHz 1% THD (BT with TV streamer)
Freq. range: 20Hz–8kHz (Mic), 20Hz–20kHz (Bluetooth)
Output power: Max 82mW @ 1kHz / 16Ω (THD=1%)

Environmental Requirements

Operating temperature: 0° – 35° C, 32° – 95° F
Relative humidity: 15% – 90% (non condensing)

FCC compliance statement

FCC ID: 2APAKBE2021

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Warning: Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Bellman & Symfon Asia Ltd.

1401 & 1408, Laifung International Centre,
33 Jie Fang Road South,
Guangzhou, China, 510115
Phone: +86 20 8318 2950
Fax: +86 20 8333 5528
E-mail: info.cn@bellman.com

Canadian compliance statement

IC: 6693A-BE2021

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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This Class B digital apparatus complies with the Canadian ICES-003.

Bluetooth® certified

This device includes a Bluetooth module which has passed the BQB certificate. The Bluetooth word mark and logos are owned by the Bluetooth SIG, Inc., and any use of such marks by Bellman & Symfon is under licenses.

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CMIIT ID: 2019DJ9374

Bellman & Symfon North America Inc.

10404 Chapel Hill Rd, Suite 102
Morrisville, NC 27560
United States
Phone: +1 (910) 239 5796
Toll-free: + (833) 235 5626
E-mail: support.bsn@bellman.com

DESIGN FOR EARS™


Manufacturer

Bellman & Symfon Group AB
Södra Långebergsgatan 30
436 32 Askim Sweden
Phone: +46 31 68 28 20
E-mail: info@bellman.com
bellman.com

Revision: BE2021-A_003MAN4.0
Date of issue: 2024-06-17
TM and © 2024 Bellman & Symfon AB.
All rights reserved.

Important safety instructions

This section contains important information about safety, handling, disposal, recycling and warranty.

Note: Read all safety information and instructions below before using the product. Keep the user manual for future use.

Warning!

- This product is intended to be operated by patient with mild to severe hearing loss.
- Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Maxi Pro BE2021, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.
- Failure to follow these safety instructions could result in fire, electric shock, or other injury or damage to the product or other property.
- Do not use or store this product near any heat sources such as naked flames, radiators, ovens or other devices that produce heat.
- Protect the cable from any potential source of damage.
- Use only power adapters and accessories specified in this manual.
- Only refer a damaged or faulty product to qualified service personnel.
- Do not dismantle the product; there is a risk of electric shock. Tampering with or dismantling the product will void warranty.
- This product is designed for indoor use only.
- Do not expose the product to moisture.
- Protect the product from shocks during storage and transport.

Care and cleaning

Note: Failure to follow these care and cleaning instructions could result in damage to the product and void the warranty.

Battery charging and replacement

Warning! The Maxi Pro is powered by a rechargeable lithium-polymer battery. Do not expose the battery to fire or to direct sunlight. The bat-

tery should only be changed by Bellman & Symfon or an authorized service center. There is a risk of explosion if the battery is not replaced correctly. The battery should be charged fully before using the product for the first time. Only charge the battery in a temperature between 0°– 35°C, 32°– 95° F.

Using the product

Operating environment: 0°– 35° C with an atmospheric pressure of 700hpa - 1060hpa. Transport and storage environment: -10°– 50° C with an atmospheric pressure of 700hpa - 1060hpa. Do not use or store the products near a heat source. If the product gets wet or is exposed to moisture, it should be replaced.

Cleaning

Use a soft, lint-free cloth. Avoid getting moisture in openings. Don't use window cleaners, household cleaners, aerosol sprays, solvents, alcohol, ammonia, or abrasives.

Service and Warranty

If the product appears to be damaged or doesn't function properly, follow the instructions in this leaflet. If the product still doesn't function as intended, contact your local dealer for information on service and warranty.

Warranty Conditions

Bellman & Symfon guarantees this product for two years from date of purchase against any defects that are due to faulty materials or workmanship. This guarantee only applies to normal conditions of use, and does not include damage resulting from accident, neglect, misuse, unauthorized dismantling, or contamination howsoever caused. This guarantee excludes incidental and consequential damage.

Model, type and classification

The information is available at the back of the product.

Product disposal and recycling information

The symbol to the left means that according to local laws and regulations your product should be disposed of separately from household waste. When this product reaches its end of life, take it to a collection point designated by local authorities.


Battery disposal information

This symbol indicates that the product contains batteries that according to local laws and regulations must be disposed of separately from household waste. Take used batteries to a collection point designated by the authorities.

Bluetooth® certified

This device includes a Bluetooth module which has passed the BQB certificate. The Bluetooth word mark and logos are owned by the Bluetooth SIG, Inc., and any use of such marks by Bellman & Symfon is under licenses.

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
Technical specifications

Intended Use

Television streamer device to be used at home in combination with the Maxi Pro assistive listening device for television volume amplification.

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Configuration options

The TV streamer can be configured with the Maxi Pro personal amplifier and the following original accessories:

- BE9228 Fiber optic cable, white
- BE9254 Single charger cable, white
- BE9273 Power adapter, white
- BE9262 Power adapter, black
- BE9123 Cable kit, white

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CMIIT ID: 2019DJ9373

Bellman & Symfon North America Inc.

5509 Business Dr, Unit B
Wilmington, NC 28405
Unites States
Phone: +1 (910) 239 5796
Toll-free: + (833) 235 5626
E-mail: support.bsn@bellman.com

DESIGN FOR EARS™

Manufacturer

Bellman & Symfon Group AB
Södra Långebergsgatan 30
436 32 Askim Sweden
Phone +46 31 68 28 20
E-mail info@bellman.com
bellman.com

Revision: BE2022-02_004MAN1.0
Date of issue: 2021-11-25
TM and © 2021 Bellman & Symfon AB.
All rights reserved.