Visit installer's guide







All you need to know about Visit.

Transmitters

5
11
15
19
23


Receivers

Flash receiver	27
Portable receiver	31
Pager receiver	35
Wrist receiver	39
Alarm clock receiver	43

Accessories

47
49
51
52
53
53
54
54

System overview





Easy access to support material

Use the product QR codes with your mobile phone or tablet to get web access to detailed product information, images and installation movies.



5 Flash receiver Alerts with bright lights



	-	
	-	
	-	
	-	

6 Alarm clock Uses sound, light & vibrations


7 Portable receiver Alerts with sound and light



8 Pager receiver Alerts with vibrations

Visit door transmitter

Buttons and connections



Technical specifications

In the box

- BE1411 Visit door transmitter
- 2 x 1.5 V AA alkaline batteries
- Velcro for wall mounting
- Screw and wall plug

Power and battery

- Battery power 2 x 1.5 V AA lithium or alkaline type batteries
- Power consumption Active < 70 mA Idle position $< 15 \,\mu$ A
- Operation time Alkaline batteries ~ 5 years Lithium batteries ~ 10 years

Dimensions and weight

- Height: 100 mm, 4.0"
- Width: 65 mm, 2.6"
- Depth: 27 mm, 1.1" Weight: 120 g, 4.2 oz. incl. batteries

Activation

- The test buttons and the int. mic.
- The electromagnetic detector
- The external microphone accessory
- The existing doorbell connected to the external trigger input

Inputs

- 3.5 mm external microphone input
- External trigger input

Environment

- For indoor use only Operating temperature 15° to 35° C, 59° to 95° F
- Relative humidity 5% to 95%, non-condensing

Frequency and coverage

- Frequency: 314.91 MHz, 433.92 MHz or 868.30 MHz, depending on region
- Coverage: 50 250 m, 55 273 yd. depending on the radio frequency and the building's characteristics


Accessories


The following accessories are available:

- BE9199 External microphone 2.5 m
- BE9200 External microphone 0.75 m

Installation – single sound source


- 1 Remove the battery pull tab to start the unit. Clean the wall with the wet wipe and remove the protective film from the Velcro.
- 2 Mount the transmitter to the left of the doorbell sound source, as close as possible. You can also use the supplied screw and plug.





Installation – intercom and electronic doorbell

- 1 Remove the protective film from the external microphone (sold separately) and attach it to the intercom speaker. Connect the external microphone to the door transmitter ext. mic. input.
- 2 Remove the battery pull tab to start the unit. Clean the wall with the wet wipe and remove the protective film from the Velcro.
- 3 Mount the transmitter to the left of the doorbell sound source, as close as possible. You can also use the supplied screw and plug.







Visit door transmitter

Installation – intercom and mechanical doorbell

If you have an intercom and a mechanical doorbell, you may need two door transmitters, i.e. one for each sound source.

- 1 Pull the battery pull tab to start each door transmitter. Clean the wall surface with the wet wipe and remove the protective film from the Velcro.
- 2 Mount one of the transmitters to the left of the doorbell and the other to the left of the intercom speaker, as close to the sound sources as possible.






Testing the connection

Using the test button

- 1 Press both test buttons simultaneously on the door transmitter. The top LED lights up in green to show that a radio signal is being transmitted.
- 2 The green Visit LED on the receiver lights up to show that the signal was received. In addition, it starts to sound, flash or vibrate depending on the receiver.

Note: If you have changed the signal pattern, it will react in accordance with the table on the following spread.




Using the doorbell or intercom

1 Press the button on the doorbell or intercom. The LED on the transmitter lights up in green to show that the sound is detected.

:::

2 The green Visit LED on the receiver lights up to show that the radio signal was received. In addition, it starts to sound, flash or vibrate depending on the receiver.

Note: If you have changed the signal pattern, it will react in accordance with the table on the following spread.




Programming the transmitter

The door transmitter is programmed from the start to recognize the majority of doorbells and intercoms. If the doorbell signal varies significantly in strength or tone, you might need to teach it to recognize the sound of your specific doorbell.

Here is how you program the transmitter:

- 1 Press and hold button 1 until the top LED starts to blink. Release it to start the recording.
- 2 Ring the doorbell at least five times. Pause 1 2 seconds between each press. The recording lasts for 5 minutes and ends automatically. You can stop it manually by pressing button 1.
- 3 When it's finished, the LED lights up in green to confirm that it was successful. If the LED is red, see Troubleshooting.
- 4 Ring the doorbell. The top LED on the transmitter lights up in green to confirm that it recognizes the sound.

You can record the sound of the intercom using button 2 by repeating the steps above.



Clearing the recorded sounds

To clear all recorded sounds, press and hold button 1 and 2 simultaneously for 5 seconds.

Note: This will clear the recorded sounds and the transmitter will return to factory settings, but the radio key settings and signal settings are kept intact.

Using electromagnetic detection

The transmitter can be set to detect electromagnetic fields generated by electric doorbells. Here is how it's done:

Open the transmitter front cover and move signal switch no. 4 to the up = on position to activate electromagnetic detection. To deactivate it, move the signal switch back to the down = off position.

Using accessories

Using an external microphone


The external microphone can be used when the sound source is located too far away from the internal microphone or when you for instance need individual notifications from the doorbell and intercom.

Connecting an accessory to the external trigger input

- Remove the front cover to access the screw terminal. There are three connection points. Use connection points 1 and 2 to connect an active switch like a relay with current. Note: The current must be 2 – 30 VDC, polarity independent or 2 – 24 VAC 5 – 150Hz.
- Use connection points 2 and 3 to connect a passive mechanical switch like a doorbell push button.









Visit door transmitter

Default signal pattern

When the door transmitter is activated, the following happens:

- 1 The LED on the transmitter lights up in green to show that it's signalling the receiver.
- 2 The Visit LED on the receiver lights up in green and it starts to sound, flash or vibrate with a certain pace, called signal pattern. The transmitter and the connected accessories determine the signal pattern. The default is as follows:

Transmitter		Receiver signal p			
Source	LED	LED	Sound	Vibration	Flash
Internal microphone	Green, top	Green light	1 x door chime, low	Slow ■□□□	Yes
 External microphone 	Green, bottom	Green blinks	2 x door chime, high	Slow ■□□□	Yes
Connected doorbell	Green, top	3 x green blinks	1 x door chime, high	Slow ■□□□	Yes

Changing the signal pattern

The transmitter controls the signal pattern. Open the transmitter front cover and move the signal switches according to the table below to change it:

On — Off —				
011	1	2	3	4

Transmit	ter	Receiver signal pattern			
Switch	Source	LED	Sound	Vibration	Flash
1 2 3 4	Int. mic. / learned signal 1 / test Ext. mic. / learned signal 2 Connected doorbell	Green light Green blink 3 x green blinks	1 x door chime, low 2 x door chime, high 1 x door chime, high	Slow Slow	Yes
1 2 3 4	Int. mic. / learned signal 1 / test Ext. mic. / learned signal 2 Connected doorbell	2 x green blinks 3 x green blinks Green blinks	2 x door chime low 1 x door chime, high 2 x door chime, high	Slow Slow	Yes
1 2 3 4	Int. mic. / learned signal 1 / test Ext. mic. / learned signal 2 Connected doorbell	3 x green blinks 2 x green blinks 3 x orange blinks	1 x door chime, high 2 x door chime, low Baby melody	Slow ■□□□ Slow ■□□□ Fast □□□□□□	Yes
1 2 3 4	Int. mic. / learned signal 1 / test Ext. mic. / learned signal 2 Connected doorbell	Green blinks Green light Orange blinks	2 x door chime, high 1 x door chime, low Baby melody	Slow Slow Fast	Yes
1 2 3 4	Int. mic. / learned signal 1 / test Ext. mic. / learned signal 2 Connected doorbell	Green light 3 x green blinks 2 x green blinks	1 x door chime, low 1 x door chime, high 2 x door chime, low	Slow ■□□□□ Slow ■□□□□ Slow ■□□□□	Yes
1 2 3 4	Int. mic. / learned signal 1 / test Ext. mic. / learned signal 2 Connected doorbell	Green light Green blinks Yellow blinks	1 x door chime low 2 x door chime, high 2 x ring signal, high	Slow ■□□□ Slow ■□□□ Medium ■□■□	Yes
1 2 3 4	Int. mic. / learned signal 1 / test Ext. mic. / learned signal 2 Connected doorbell	3 x green blinks Green light Red + orange light	1 x door chime, high 1 x door chime, low Emergency siren	Slow ■□□□ Slow ■□□□ Long ■□□□	Yes
1 2 3 4	Int. mic. / learned signal 1 / test Ext. mic. / learned signal 2 Connected doorbell	Green blinks 2 x green blinks Red blinks	2 x door chime, high 2 x door chime, low Fire horn	Slow Slow	Yes


Changing the radio key

If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference, you need to change the radio key on all units. The radio key switches are located under the transmitter cover.

Here is how you change the radio key:

- 1 Open the transmitter front cover and move any radio key switch to the up = on position to change the radio key. (By default, all radio key switches are positioned down = off.)
- 2 Press and hold the test/function button on the receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- 3 Press both test buttons simultaneously on the transmitter within 30 s to send the new radio key.
- 4 All Visit LEDs on the receiver blink 5 times to show that the radio key has been changed. It then returns to normal mode.

Please note: All Visit units must be set to the same radio key in order to operate as a group.



Troubleshooting

If	Try this
The LEDs blink in orange every minute	 Replace the batteries
The transmitter LEDs blinks in orange every second	 There are other comp them off or turn then Disconnect the extern
The transmitter LED lights up when I press the doorbell or intercom – but the receiver is not activated	 Check the the transm Move the receiver clc Check that the door t For more information
The transmitter LED doesn't light up when I ring the doorbell or intercom	 Ring the doorbell wh the sound source. The Program the transmit If the signal varies a lo If the door transmittee simultaneously for 5 statements
The receiver is activated for no apparent reason	 There is probably and system. Change the r






s. Only use 1.5 V AA (LR6) lithium or alkaline batteries.


- peting sound sources around the door transmitter. Switch m down.
- rnal microphone accessory to make sure it is not faulty.
- nitter batteries and the receiver batteries and connections. oser to the transmitter to make sure it's within radio range. transmitter and the receiver are set to the same radio key. n, see Changing the radio key.
- hile moving the transmitter closer and further away from ne ideal distance is less than 3 cm.
- itter to recognize the doorbell sound. See Programming. lot in strength or tone, change to electromagnetic detection. ter is still not activated, press and hold button 1 and 2 s to clear the recorded sounds and repeat the steps above.

other Visit system installed nearby that triggers your radio key on all units, see Changing the radio key.

Visit push button transmitter

Buttons and connections





Technical specifications

In the box

- BE1420 Visit push button transmitter with pre-mounted alkaline battery
- Lanyard with safety clasp
- Adhesive tape, screws and plugs

Power and battery

- Battery type 1 x 6 V PX28A alkaline or 1 x 6 V PX28L lithium
- Power consumption Active < 35 mA Idle position < 0.05 μA
- Operation time Alkaline battery ~ 2 years Lithium battery ~ 5 years

Dimensions and weight

- Height: 66 mm, 2.6"
- Width: 48 mm, 1.9"
- Depth: 23 mm, 0.9"
- Weight: 50 g, 1.8 oz. incl. battery

Activation

Via the front button

Maintenance and cleaning

- Maintenance free Clean with a dry cloth
- Do not use household cleaners, aerosol sprays, solvents, alcohol, ammonia or abrasives

Frequency and coverage

- Frequency 314.91 MHz, 433.92 MHz or 868.30 MHz, depending on the region
- Coverage 50 - 250 m, 55 - 273 yd. depending on the radio frequency and the characteristics of the building

Environment

- For indoor use and outdoor use in a protected location. Will not withstand water or rain.
- Operating temperature 15° to 35° C, 59° to 95° F
- Relative humidity 5% to 95%, non-condensing

Using it as a caller button

The transmitter can be worn around your neck and be used as a wireless caller button.

Here is how you use it:


- 1 Attach the lanyard to the transmitter.
- 2 Hang the transmitter around your neck.

Using it as a doorbell

The transmitter can also be used as a wireless doorbell. Here is how you set it up:


Mounting with screws

- 1 Carefully remove the transmitter front cover and make two holes on the markings.
- 2 Fix the unit to the wall using the supplied screws and put the cover back on.




Mounting with adhesive tape


- 1 Attach the adhesive tape to the back of the transmitter. Clean the wall with the wet wipe and remove the protective film from the tape.
- 2 Mount the unit in a weather protected area by the front door.












Visit push button transmitter

Testing the connection

Using the front button

- 1 Press the front button on the transmitter. The LED lights up in green to show that a radio signal is being transmitted.
- 2 The green Visit LED on the receiver lights up to show that the signal was received. In addition, it starts to sound, flash or vibrate with a certain pace, called signal pattern.

The transmitter determines the signal pattern and the default is as follows:



Or

Off

SIGNAL

Default signal pattern

Transmitter	Receiver signal pattern						
LED	LED	Sound		Vibrati	on		Flash
 Green light 	Green light	1 x door chime, low		Slow			Yes

Changing the signal pattern

Changing the signal pattern is easy. Just open the transmitter front cover and move signal switches no. 8, 9 and 0 according to the table below:

Transmitter **Receiver signal pattern**

Switch	LED	Sound	Vibration	Flash
890	Green light	1 x door chime, low	Slow ■□□□	Yes
8 9 O	2 x green blinks	2 x door chime, low	Slow ■□□□	Yes
890	3 x green blinks	1 x door chime, high	Slow ■□□□	Yes
890	Green blinks	2 x door chime, high	Slow ■□□□	Yes
890	Orange light	Baby melody	Fast IDIDID	Yes
890	Orange blinks	Baby melody	Fast Intoin	Yes
890	Yellow light	1 x ring signal, low	Medium ∎□∎□	Yes
890	Yellow blinks	2 x ring signal, high	Medium ∎□∎□	Yes


Changing the radio key

If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference you need to change the radio key on all units. The radio key switches are located under the transmitter cover.

Here is how you change the radio key:

- 1 Remove the transmitter front cover and move any radio key switch to the up = on position to change the radio key. (By default, all radio key switches are positioned down = off.)
- 2 Press and hold the test/function button on the receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- **3** Press the front button on the transmitter within 30 seconds to send the new radio key.
- 4 All Visit LEDs on the receiver blink 5 times to show that the radio key has been changed. It then returns to normal mode.

Note: All Visit units must be set to the same radio key in order to operate as a group.






Troubleshooting

If	Try this
The transmitter LED lights up in yellow when I press the button	 The battery is nearly of PX28L type battery.
The transmitter LED doesn't light up when I press the button	Check that the batterReplace the battery w
The transmitter LED lights up in green but the receiver is not activated	 Check the receiver ba Move the receiver clo Check that the units a
The receiver is activated for no apparent reason	 There is probably and system. Change the r








depleted. Replace it with an alkaline PX28A or a lithium

- ery is positioned correctly. with an alkaline PX28A or a lithium PX28L type battery.
- atteries and connections.
- loser to the transmitter to make sure it's within radio range. are set to the same radio key, see Changing the radio key.
- nother Visit system installed nearby that triggers your radio key on all units, see Changing the radio key.

Visit telephone transmitter

Buttons and connections



Technical specifications

In the box

- BE1431 Visit telephone transmitter
- 2 x 1.5 V AA alkaline batteries
- Telephone cord and adapter
- Screw and wall plug

Power and battery

- Battery power
 2 x 1.5 V AA lithium or alkaline
 type batteries
- Power consumption Active < 70 mA Idle position < 15 μA
- Operation time Alkaline batteries ~ 5 years Lithium batteries ~ 10 years

Dimensions and weight

- Height: 100 mm, 4.0"
- Width: 65 mm, 2.6"
- Depth: 27 mm, 1.1"
- Weight: 120 g, 4.2 oz. incl. batteries

Activation

- The test buttons
- A landline telephone
- A smartphone or tablet via the mobile phone sensor
- A contact mat or magnetic switch
- A doorbell connected to the ext trig

Environment

For indoor use only

Accessories

- BE9251 Mobile phone sensor
- BE9023 Magnetic switch
- BE9024 Contact mat
- BE9253 Ext. trig. cable, 3.5 mm

Inputs

- RJ11 analogue telephone input
- 2.5 mm external trigger input
- 3.5 mm external trigger input


Frequency and coverage


- Frequency: 314.91 MHz, 433.92 MHz or 868.30 MHz, depending on region
- Coverage: 50 250 m, 55 273 yd. depending on the radio frequency and the building's characteristics

Setting up the transmitter

- Remove the battery pull tab to start the unit.
 Clean the wall with the wet wipe and remove the protective film from the Velcro.
- 2 Mount the transmitter on the wall. You can also use the supplied screw and plug.
- **3** Connect the telephone adapter as shown below.

Note: The appearance of the adapter may differ with territory.






Testing the connection

Using the test button

- 1 Press both test buttons simultaneously on the telephone transmitter. The top LED lights up in green to show that a radio signal is being transmitted.
- 2 The yellow Visit LED on the receiver lights up to show that the signal was received. In addition, it starts to sound, flash or vibrate depending on the receiver.

Note: If you have changed the signal pattern, it will react in accordance with the table on the following spread.








film from the Velcro.

Using a mobile phone

- Use for instance a mobile phone to call the landline telephone. The top LED on the transmitter lights up in green to show that an incoming call is detected.
- 2 The yellow Visit LED on the receiver lights up to show that the radio signal was received. In addition, it starts to sound, flash or vibrate depending on the receiver.

Note: If you have changed the signal pattern, it will react in accordance with the table on the following spread.



Visit telephone transmitter

Default signal pattern

When the telephone transmitter is activated by an incoming call or a triggered accessory, the following happens:

- 1 The LED on the transmitter lights up to show that it's signalling the receiver.
- 2 The Visit LED on the receiver lights up and it starts to sound, flash or vibrate with a certain pace, called signal pattern. The transmitter and the connected accessories determine the signal pattern. The default is as follows:

Transmitter		Receiver signal pattern			
Source	LED	LED	Sound	Vibration	Flash
Landline phone	Green, top	Yellow light	1 x ring signal, low	Medium ∎□∎□	Yes
 Mobile phone sensor 	Green, top	Yellow blinks	2 x ring signal, high	Medium ∎□∎□	Yes
 Other accessory 	Green, bottom	Green light	1 x door chime, low	Slow ■□□□	Yes

Changing the signal pattern

The transmitter controls the signal pattern. Open the transmitter front cover and move the signal switches according to the table below to change it:

On – Off –	- 1	
011	123	4
	SIGNA	

Transmitter		Receiver signal pattern			
Switch	Source	LED	Sound	Vibration	Flash
	Landline phone / test button Mobile phone sensor	Yellow light Yellow blinks	1 x ring signal, low 2 x ring signal, high	Medium ∎□∎□ Medium ∎□∎□	Yes Yes
1234	Other accessory	Green light	1 x door chime, low	Slow ■□□□	Yes
	Landline phone / test button	Yellow light	1 x ring signal, low	Medium ∎□∎□	Yes
	Mobile phone sensor	Yellow blinks	2 x ring signal, high	Medium ∎□∎□	Yes
1 2 3 4	Other accessory	2 x green blinks	2 x door chime, low	Slow ■□□□	Yes
	Landline phone / test button	Yellow light	1 x ring signal, low	Medium ∎□∎□	Yes
	Mobile phone sensor	Yellow blinks	2 x ring signal, high	Medium ∎□∎□	Yes
1 2 3 4	Other accessory	3 x yellow blinks	1 x ring signal, high	Medium ∎□∎□	Yes
	Landline phone / test button	Yellow light	1 x ring signal, low	Medium ∎□∎□	Yes
	Mobile phone sensor	Yellow blinks	2 x ring signal, high	Medium ∎□∎□	Yes
1 2 3 4	Other accessory	2 x orange blinks	Baby melody	Fast Internet	Yes
	Landline phone / test button	2 x yellow blinks	2 x ring signal, low	Fast Information	Yes
	Mobile phone sensor	Yellow light	1 x ring signal, low	Medium ∎□∎□	Yes
1 2 3 4	Other accessory	3 x orange blinks	Baby melody	Fast Internet	Yes
	Landline phone / test button	2 x yellow blinks	2 x ring signal, low	Medium ∎□∎□	Yes
	Mobile phone sensor	Orange blinks	Baby melody	Fast Internet	Yes
1234	Other accessory	Green blinks	2 x door chime, high	Slow ■□□□	Yes
	Landline phone / test button	Orange blinks	Baby melody	Medium ∎□∎□	Yes
	Mobile phone sensor	3 x yellow blinks	1 x ring signal, high	Medium ∎□∎□	Yes
1 2 3 4	Other accessory	2 x green blinks	2 x door chime, low	Slow ■□□□	Yes
	Landline phone / test button	3 x yellow blinks	1 x ring signal, high	Medium ∎□∎□	Yes
	Mobile phone sensor	2 x yellow blinks	2 x ring signal, low	Medium ∎□∎□	Yes
1 2 3 4	Other accessory	Green blinks	2 x door chime, high	Slow	Yes


Changing the radio key

If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference you need to change the radio key on all units. The radio key switches are located under the transmitter cover.

Here is how you change the radio key:

- 1 Open the transmitter front cover and move any radio key switch to the up = on position to change the radio key. (By default, all radio key switches are positioned down = off.)
- 2 Press and hold the test/function button on the receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- 3 Press both test buttons simultaneously on the transmitter within 30 seconds to send the new radio key.
- 4 All Visit LEDs on the receiver blink 5 times to show that the radio key has been changed. It then returns to normal mode.

Note: All Visit units must be set to the same radio key in order to operate as a group.




Troubleshooting

If	Try this
The LEDs blink in orange every minute	Replace the batteries.
The transmitter LED lights up in green but the receiver doesn't respond	 Check the the transm Move the receiver clo Check that the units a
The transmitter LED doesn't light up when the phone rings or when an accessory is triggered	 Press the test buttons connections. If the LE 1.5 V AA (LR6) lithium
The transmitter LED doesn't light up when I press the test buttons	 Replace the batteries. LED still doesn't light
The receiver is activated for no apparent reason	 There is probably and system. Change the rate









s. Only use 1.5 V AA (LR6) lithium or alkaline batteries.

nitter batteries and the receiver batteries and connections. oser to the transmitter to make sure it's within radio range. are set to the same radio key, see Changing the radio key.


is on the transmitter. If the LED lights up in green, check all ED doesn't light up in green, replace the batteries. Only use n or alkaline batteries.

s. Only use 1.5 V AA (LR6) lithium or alkaline batteries. If the t up, contact your retailer for service information.

other Visit system installed nearby that triggers your radio key on all units, see Changing the radio key.

Visit baby monitor

Buttons and connections



Always make sure the baby monitor is out of the child's reach. Never place the baby monitor in the child's crib or playpen.

Technical specifications

In the box

- BE1491 Visit baby monitor
- 2 x 1.5 V AA alkaline batteries
- Pre-mounted table stand
- Screw and wall plug

Power and battery

- Battery power 2 x 1.5 V AA lithium or alkaline type batteries
- Power consumption Active < 70 mA Idle position $< 400 \,\mu$ A
- Operation time Alkaline batteries ~ 6 months Lithium batteries ~ 1 year

Dimensions and weight

- Height: 100 mm, 4.0"
- Width: 65 mm, 2.6"
- Depth: 27 mm, 1.1"
- Weight: 120 g, 4.2 oz. incl. batteries

Activation

- The internal mic. and test buttons
- The contact mat accessory

Settings

- Sensitivity 65 dB, 75 dB, 85 dB
- Delay 30 sec, 10 sec, 1 sec

Environment

BE9024 contact mat

- For indoor use only Operating temperature 15° to 35° C, 59° to 95° F
- Relative humidity 5% to 95%, non-condensing

Frequency and coverage


- Frequency: 314.91 MHz, 433.92 MHz or 868.30 MHz, depending on region
- Coverage: 50 250 m, 55 273 yd. depending on the radio frequency and the building's characteristics

Accessories

BE9024 Contact mat Alerts you if your baby leaves the bed

Setting up the baby monitor

- 1 Remove the battery pull tab to start the unit. (You can press and hold both test buttons for 3 seconds to turn it on/off.)
- 2 Place the baby monitor on the bedside table or mount it on the wall using the supplied screw and plug. The recommended distance is 0.5 – 2 m, always out of reach from the child.




Testing the connection

Using the test button

- 1 Press both test buttons simultaneously on the baby monitor. The top LED lights up in green to show that a radio signal is being transmitted.
- 2 The orange Visit LED on the receiver lights up to show that the signal was received. In addition, it starts to sound, flash or vibrate depending on the receiver.

Note: If you have changed the signal pattern, it will react in accordance with the table on the following spread.






Using your voice

- 1 Stand by the baby monitor and make some noise. The top LED lights up in green to show that the sound was detected.
- 2 The orange Visit LED on the receiver lights up to show that the radio signal was received. In addition, it starts to sound, flash or vibrate depending on the receiver.

Note: If you have changed the signal pattern, it will react in accordance with the table on the following spread.



Visit baby monitor

Adjusting the settings

The buttons for sensitivity and delay are located on the left side of the baby monitor. When you press the button, the corresponding LED colour shows the current setting. Then press repeatedly to adjust the setting.


- If the baby monitor is not activated when the baby cries increase the sensitivity.
- If the baby monitor is activated too easily reduce the sensitivity.
- If the baby monitor is activated too early or too late adjust the delay.

Sensitivity settings

Delay settings



30 s Red Ö 10 s Orange Green 1 s



Default signal pattern

When the baby monitor is activated by the baby's voice or the contact mat, the following happens:

- 1 The LED on the baby monitor lights up in green to show that it is signalling the receiver.
- 2 The Visit LED on the receiver lights up in orange and it starts to sound, flash or vibrate with a certain pace, called signal pattern. The baby monitor and the contact mat accessory determine the signal pattern. The default is as follows:

Receiver signal pattern Baby monitor Source LED LED Sound Vibration

 Baby voice 	Green, top	Orange light	Baby melody	Fast Interne	Yes
Contact mat	Green, bottom	Green light	1 x door chime, low	Slow ∎□□□	Yes

Changing the signal pattern

Contact mat 4

If you have more than one child, you can set a unique Visit LED pattern for each baby monitor. Just open the front cover and move any signal switch to the up = on position and make sure the other switches are down = off.

On — Off —	
	1 2 3 4

Yes

Fast Internet

Flash

SIGNAL **Baby monitor Receiver signal pattern** Switch Source LED Sound Vibration Flash Baby monitor 1 Orange light Baby melody Fast Information Yes Contact mat 1 Green light Baby melody Fast IDIDIDID Yes Baby monitor 2 2 x orange blinks Baby melody Fast IDIDIDID Yes Contact mat 2 2 x green blinks Yes Baby melody Fast Information Baby monitor 3 3 x orange blinks Baby melody Yes Fast Information 1234 Contact mat 3 3 x orange blinks Baby melody Fast IDIDIDID Yes Orange blinks Baby monitor 4 Baby melody Yes Fast Information

Orange blinks

Baby melody


Changing the radio key

If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference, you need to change the radio key on all units. The radio key switches are located under the monitor cover.

Here is how you change the radio key:

- 1 Carefully remove the table stand and open the front cover on the baby monitor. Move any radio key switch to the up = on position to change the radio key.
- 2 Press and hold the test/function button on the receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- 3 Press both test buttons simultaneously on the baby monitor within 30 s to send the new radio key.
- 4 All Visit LEDs on the receiver blink 5 times to show that the radio key has been changed. It then returns to normal mode.

Note: All Visit units must be set to the same radio key in order to operate as a group.



Troubleshooting

lf	Try this
The LEDs blink in orange every minute	 Replace the batteries.
The baby monitor LED lights up in green but the receiver is not activated	 Check the baby monit Move the receiver close Check that the units a
The baby monitor LED doesn't light up even though the baby is crying	 Move the baby monit sensitivity, see Adjust Note: Always m Never place the
The baby monitor is activated too easily	 Reduce the sensitivity Adjusting the setting
The baby monitor is activated too early	Increase the delay. For
The baby monitor is activated too late	 Reduce the delay, see
The receiver is activated for no apparent reason	 There is probably ano system. Change the ratio

m





On


Off

. Only use 1.5 V AA (LR6) lithium or alkaline batteries.

- itor batteries and the receiver batteries and connections. oser to the baby monitor to make sure it's within range. are set to the same radio key, see Changing the radio key.
- itor closer to the baby or increase the microphone sting the settings.
- nake sure the baby monitor is out of the child's reach. e baby monitor in the child's crib or playpen.
- y or move the baby monitor further away from the crib, see gs.
- or more information, see Adjusting the settings.
- e Adjusting the settings.
- other Visit system installed nearby that triggers your radio key on all units, see Changing the radio key.

Visit smoke alarm

Buttons and controls



Technical specifications

In the box

- BE1480 Visit smoke alarm
- 1 x 9 V alkaline or lithium battery
- Mounting plate
- Screws and plugs

Power and battery

- Battery power 9V Duracell MN1604, Energizer 522 or 9V Ultralife U9VL-J (lithium)
- Power consumption Active < 40 mA Idle position $< 10 \,\mu$ A
- Operation time Alkaline battery ~ 3 years Lithium battery ~ 6 years

Dimensions and weight

- Height: 100 mm, 3.9"
- Width: 100 mm, 3.9"
- Depth: 35 mm, 1.4"
- Weight: 110 g, 3.9 oz. incl. battery

Sensor type

Optothermal sensor with an audible alarm of > 85 dB(A) @ 3m

Activation

- The test button
- The built-in smoke detector
- The built-in temperature sensor, if the temperature exceeds ~ 57°C.

Environment

- For indoor use only Operating temperature 0° to 38° C, 32° to 100° F
- Relative humidity 15% to 95%, non condensing

Frequency and coverage


- Frequency: 314.91 MHz, 433.92 MHz or 868.30 MHz, depending on region
- Coverage: 50 250 m, 55 273 yd. depending on the radio frequency and the building's characteristics

Regulatory requirements

BE1480 complies with the smoke alarm standard EN 14604:2005

Setting up the smoke alarm

- 1 Remove the mounting plate and connect the battery to the battery snaps to start the unit. Wait for about 10 seconds while the smoke alarm carries out a self-test. It is finished when the test button blinks once.
- 2 Fix the mounting plate to the ceiling using the supplied screws and plugs. Mount it at least 50 cm, (19.7") from walls and other obstructions, see Fitting the smoke alarm.
- 3 Attach the smoke alarm to the mounting plate by turning it clockwise. Extend the radio antenna so that it points down.




Testing the connection

Using the test button

- 1 Press and hold the test button on the smoke alarm for more than one second. The smoke alarm will beep and transmit a radio signal to the receiver.
- 2 The red Visit LED on the receiver lights up to show that the signal was received. In addition it starts to sound, flash or vibrate depending on the receiver.

Note: Once you release the smoke alarm test button, the beep will time out in a couple of seconds.






Using smoke detector test aerosol

- 1 Spray some of the test material into the chamber and wait 5 - 10 seconds for the smoke alarm to beep and transmit a radio signal to the receiver.
- 2 The red Visit LED on the receiver lights up to show that the signal was received. In addition it starts to sound, flash or vibrate depending on the receiver.

Note: The smoke alarm will beep and transmit the signal as long as there is test aerosol inside the chamber.



Visit smoke alarm

Default signal pattern

The smoke alarm LED blinks in red once per minute to show that it is working correctly. Depending on the alarm, the signal patterns are as follows:

Smoke alarm			Receiver signal pattern				
Alarm type	LED	Sound	LED	Sound	Vibrati	ion	Flash
 Fire detected 	Red blinks	Fire alarm	Red blinks	Fire alarm	Long		Yes
Low battery	2 x red blinks	1 beep every min	Red blink every 5 s	None	None		None
 Flat battery 	Red blinks	1 beep every s	Red blink every 5 s	1 x fire alarm	Slow		None

Using broadcast mode

If you want the smoke alarm signal to be transmitted to *all* Visit receivers within radio range, you can activate broadcast mode. This will override the radio key settings.

Here is how you activate broadcast mode:

Move radio switch 7 on the back of the smoke alarm to the up = on position to activate broadcast mode.

Note: Activation with the test button and battery warnings will only be transmitted to units with the same radio key.

Using toast mode

To avoid false alarm when you are for instance cooking, you can temporarily reduce the smoke alarm sensitivity.

Here is how you activate toast mode:

- Press the test button briefly. The smoke alarm will beep and blink twice in yellow. The LED will continue to blink in red while it's in toast mode.
- Press the test button again to deactivate toast mode. The smoke alarm will beep and blink three times in yellow.

Note: The toast mode times out in 20 minutes.

Fitting the smoke alarm

Fit the smoke alarm in the center of the ceiling outside the bedrooms, at least 50 cm from any wall. If the bedrooms are in different areas of the house, separate smoke alarms are recommended. In multi-storey properties, install at least one smoke alarm on each floor.

Avoid kitchens, fireplaces or garages, as cooking fumes and car exhaust may cause false alarms. The smoke alarm should not be installed in damp spaces, close to fans, etc. or in agricultural buildings.

Testing and maintenance

Test the smoke alarm regularly, preferably each week, e.g. during cleaning, but at least once per month. Always test it immediately after any holidays or other extended periods of absence. Clean it with a damp cloth. After changing battery, vacuum clean with a soft brush. Do not paint over the smoke alarm.

Minimum, O = Additional


Changing the radio key

If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference you need to change the radio key on all units. The radio key switches are located on the back of the smoke alarm.

Here is how you change the radio key:

- 1 Remove the mounting plate and move any of the radio key switches 1 6 to the up = on position to change the radio key.
- 2 Press and hold the test/function button on the receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- 3 Press the test button on the smoke alarm for more than one second to send the new radio key.
- 4 All Visit LEDs on the receiver blink 5 times to show that the radio key has been changed. It then returns to normal mode.

Note: All Visit units must be set to the same radio key in order to operate as a group. If broadcast mode is activated, all Visit receivers will respond regardless of the radio key settings.




Troubleshooting

Most problems with the smoke alarm can be solved quickly by following the advice below.

lf	Try this
Nothing happens when I press the test button	 Replace the smoke al Use a Duracell MN16
The smoke alarm beeps when I press the test button, but the receiver is not responding	 Check that the smoke Check the smoke alar Move the receiver close Check that the units a
The receiver is activated for no apparent reason	 Replace the battery. If the problem persist triggers yours. Change
The receiver beeps and chirps for no apparent reason	 The smoke alarm is d







alarm battery. 604, Energizer 522 or Ultralife U9VL-J type battery.

ke alarm antenna is straight and points to the floor.

arm battery and the receiver batteries and connections.

loser to the smoke alarm to make sure it's within radio range.

are set to the same radio key, see Changing the radio key


Use a Duracell MN1604, Energizer 522 or Ultralife U9VL-J. sts, there is probably another Visit system nearby that ge the radio key on all units, see Changing the radio key.

defective and needs to be sent for repair.

Visit flash receiver

Buttons and controls

Visit LEDs Flash Mute / test button 8 🛽 🖊 🛦 • • • • 🗲 🖉 Bellm BE1442 batteries Cable holder Power LED



1 - Power supply, 2 - BE9251 mobile phone sensor 3 - BE1270 bed shaker, 4 - BE9105 telephone cable

Technical specifications

In the box

- BE1441 Flash receiver or BE1442 Flash receiver w. battery backup
- Power supply
- 4 x 1.2 V AAA NiMH batteries (BE1442 model only)

Power and battery

- Mains power 7.5 V DC / 1500 mA External power supply unit
- Power consumption Active: 1250 mA, idle position: 10 mA
- Backup batteries (BE1442 model only) 4 x 1.2 V AAA NiMH rechargeable batteries
- Backup battery operating time ~ 48 h when fully charged

Accessories

Dimensions and weight

Diameter BE1441: 70 mm, 2.7"

The Visit LEDs normally indicates

Orange LED, pacifier symbol

Green LED, door symbol

Red LED, fire symbol

The baby monitor is activated

The door transmitter is activated

The phone transmitter is activated

Yellow LED, telephone symbol

The smoke alarm is activated

BE1441: 140 mm, 5.5"

BE1442: 155 mm, 6.1"

BE1442: 78 mm, 3.1"

310 g, 10.9 oz.

Height

Weight

Visit LEDs

the following:

- BE9075 Wall bracket
- BE1270 Bed shaker
- BE9251 Mobile phone sensor
- BE9105 Telephone cable

Frequency and coverage


- Radio frequency 314.91 MHz, 433.92 MHz or 868.30 MHz, depending on the region
- Coverage 50 - 250 m, 55 - 273 yd. depending on the radio frequency and the characteristics of the building

Output

Built-in ~30 Candela Xenon light Warning! Flashes can cause epileptic attacks

Getting started

- 1 Connect the power supply to the receiver and the mains outlet. Pull the battery tab on the bottom (BE1442 only.) Place the receiver on a level surface or mount it on the wall using the wall bracket accessory (see separate instructions).
- 2 To test the radio link you need a Visit transmitter. Press the test button/s on the transmitter.
- 3 The receiver lights up a Visit LED and starts to flash. If a bed shaker is connected, it will vibrate. A short press on the mute/test button repeats the last indication. If nothing happens, see Troubleshooting.



Default signal pattern

When a transmitter is activated, the flash receiver lights up a LED, flashes and the bed shaker starts to vibrate with a certain pace. This is called signal pattern. The transmitters determine the pattern, and the default is as follows:

Transmitter Flash Activated source Visit L Door transmitter / push button transmitter Green Telephone transmitter / connected telephone Yellow

- Baby monitor Orang
- Smoke alarm Red

Changing the signal pattern

The signal pattern can only be changed on the transmitters. See Changing the signal pattern for the relevant transmitter.

Power LED indications

When the flash receiver is connected to mains power, the power LED lights up in green. The BE1442 model is also equipped with battery backup and the power LED indicates the following:

Status

Power LED

- Green light
- Green blinks
- Red light
- Red blinks
- The flash receiver is connected to mains power. The backup batteries are detected. The flash receiver is connected to mains power. No backup batteries are detected. The flash receiver is running on battery backup.
- The backup batteries are nearly depleted.

27



receiver		Bed shaker
LED	Flash	Vibration
า	Yes	Slow ■□□□
N	Yes	Medium ∎□∎□
ge	Yes	Fast IDIDIDID
	Yes	Long


Visit flash receiver

Changing the radio key

If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference, you need to change the radio key on all units. The radio key switches are located on the transmitters.

Here is how you change the radio key:


1 Open the transmitter cover and move any radio key switch to the up = on position to change the radio key. See Changing the radio key for the relevant transmitter.



- 3 Press the test button/s on the transmitter within 30 seconds to send the new radio key.
- 4 All Visit LEDs on the receiver blink 5 times to show that the radio key has been changed. It then returns to normal mode.

Off

Please note: All Visit units must be set to the same radio key in order to operate as a group.



Accessories

The flash receiver can be complemented with the following accessories:

BE1270 Bed shaker

Wakes you with vibrations under the pillow or mattress.

BE9250 Mobile phone sensor

Place it on the mobile phone or tablet, and the flash receiver will alert you of incoming calls and messages.

BE9105 Telephone cable

Use it to connect the receiver to your landline telephone and be alerted with flashes when the phone rings.

BE9075 Wall bracket

Directing the flash

The flash receiver features a rotating top that makes it easy to direct the light. Point it for example towards a wall if you feel that the flash is too intense. A silicone slip-on top is also available in a variety of colours (art. no. BE9164-BE9167).

Advanced programming

By using advanced programming, you can customize the signal pattern from a specific transmitter and event, displaying the LED colour and vibration pattern of your choice. The advanced programming overrides the radio key and pairs the units via the serial number. Please note that smoke alarms cannot be programmed for safety reasons.

Note: The transmitter must be activated as it is intended to be used in the system to generate the right signal. This means that you can't always use the transmitter test button (see **Default signal pattern** for the relevant transmitter).

Here is how you program the receiver:

- 1 Press and hold the mute/test button on the receiver. The green and yellow Visit LEDs will start to blink alternately. While still holding down the button, activate the desired transmitter as intended. Release the button.
- 2 Scroll through the different Visit LED options by pressing the mute/test button on the receiver. Select the desired Visit LED colour by holding down the mute/test button until the power LED goes out and lights up again.
- 3 Scroll through the different vibration options by pressing the test button on the receiver (bed shaker required). Select the desired vibration pattern by holding down the mute/test button until the power LED goes out and lights up again.
- 4 The receiver will now show the new Visit LED colour and vibration pattern. Press the mute/test button briefly to end the demonstration. After a short while, it will return to normal mode.

Deleting the advanced programming

Follow the procedure below to delete the advanced programming.

- 1 Hold down the mute/test button on the receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- 2 Press the mute/test button on the receiver 3 times in quick succession.
- 3 All Visit LEDs will light up for ~2 seconds to show that it has been deleted.

Troubleshooting

Most problems with the flash receiver can be solved quickly by following the advice below.

If	Try this
The receiver seems to be turned off	 Check that the power sup Charge the backup batter
The power LED blinks in red	 The backup batteries are Connect the power suppl
The power LED blinks in green	 The receiver detects no b
The receiver does not respond when a transmitter is activated, but works when I use the test button	 Check the transmitter bat Move the receiver closer t Check that the receiver is system, see Changing the
The receiver is activated for no apparent reason	 There is probably another Change the radio key on a
The flash is too bright	 Redirect the light by rotat



BE1441

pply is connected correctly. ries for at least 24 hours (BE1442 only).

nearly depleted and the power supply is disconnected. ly and charge the batteries for at least 24 hours.

backup batteries. Pull the battery tab, see Getting started.

tteries and connections.


to the transmitter to make sure it's within radio range. s set to the same radio key as the other units in the Visit e radio key.

er Visit system installed nearby that triggers your system. all units, see Changing the radio key.

ting the top or use a silicone slip-on top to dim the light.

Visit portable receiver

Buttons and controls



Technical specifications

In the box

- BE1450 Portable receiver
- 4 x 1.5V LR14 batteries
- Wall bracket
- Screws and wall plugs

Power and battery

- Mains power 7.5 V DC / 1000 mA Optional power supply unit Europe: BE9201, UK: BE9202
- Battery power 4 x 1.5 V LR14 alkaline batteries
- Operating time 2 – 3 years with alkaline batteries
- Power consumption Active: 1000 mA, idle position: 0.1 mA

Dimensions and weight

- Height: 165 mm, 6.5" Width: 130 mm, 5.1"
- Depth: 36 mm, 1.5"
- Weight: 590 g, 20.8 oz. incl. batteries

Visit LEDs

The Visit LEDs normally indicates the following:

- Orange LED, pacifier symbol The baby monitor is activated
- Green LED, door symbol The door transmitter is activated
- Yellow LED, telephone symbol The phone transmitter is activated
- Red LED, fire symbol The smoke alarm is activated

Output

- Adjustable sound signal Max 93 dBA @ 1 m, frequency range: 500 – 1000 Hz
- Bed shaker outlet: 2.0 4.0 VDC or Speaker outlet: $10 \text{ k}\Omega$, 0 - 4 V

Frequency and coverage


- Radio frequency 314.91 MHz, 433.92 MHz or 868.30 MHz, depending on the region
- Coverage 50 - 250 m, 55 - 273 yd. depending on the radio frequency and the characteristics of the building

Accessories

- BE1270 Bed shaker
- BE9201/BE9202 Power supply unit

Getting started

- 1 Slide open the battery cover, fit the batteries and close the cover again. Place the receiver on a level surface or mount it on the wall using the wall bracket.
- 2 To test the radio link you need a Visit transmitter. Press the test button/s on the transmitter.
- 3 The receiver lights up a Visit LED and starts to flash and sound. If a bed shaker is connected, it will vibrate. A short press on the test button repeats the last indication. If nothing happens, see Troubleshooting.



Default signal pattern

When a transmitter is activated, the receiver lights up a LED, sounds, flashes and the bed shaker starts to vibrate with a certain pace. This is called signal pattern. The transmitters determine the pattern, and the default is as follows:

Transmitter	Portable re	Portable receiver		
Activated source	Visit LED	Sound	Flash	Vibration
 Door transmitter / push button transmitter 	Green	Door chime	Yes	Slow ■□□□
 Telephone transmitter 	Yellow	Ring signal	Yes	Medium ∎□■□
 Baby monitor 	Orange	Baby melody	Yes	Fast IDIDIDID
Smoke alarm	Red	Fire horn	Yes	Long

Changing the signal pattern

The signal pattern can only be changed on the transmitters. See Changing the signal pattern for the relevant transmitter.

Adjusting the volume and flash


Adjust the volume to your liking using the red volume dial on the top of the receiver. It goes from 0 to 93 dBA @ 1 m with a main frequency range of 500 – 1000 Hz. Use the flash signal switch on the back of the receiver to turn the flash off/on.

Replacing batteries

If the power LED is yellow when the receiver is activated, the batteries are nearly depleted. Here is how you replace them: Slide open the battery cover. Replace the old batteries with four new 1.5 V LR14 alkaline batteries, see the battery

compartment for correct positioning.






Visit portable receiver

Changing the radio key

If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference, you need to change the radio key on all units. The radio key switches are located on the *transmitters*.


Here is how you change the radio key:

1 Open the transmitter cover and move any radio key switch to the up = on position to change the radio key. See Changing the radio key for the relevant transmitter.



- 3 Press the test button/s on the transmitter within 30 seconds to send the new radio key.
- 4 All Visit LEDs on the receiver blink 5 times to show that the radio key has been changed. It then returns to normal mode.

Please note: All Visit units must be set to the same radio key in order to operate as a group.



Accessories

The portable receiver can be complemented with the following accessories:


BE1270 Bed shaker

Wakes you with vibrations if anything happens while you are asleep. Connect it to the receiver and slide it under your pillow or mattress.

BE9201 EU / BE9202 UK power supply

If your receiver has a fixed place, you can connect it to mains power and not having to worry about batteries.

Connecting a speaker



Off

On the receiver model BE1452, the bed shaker outlet is replaced with a 3.5 mm audio output that can drive an external speaker. The speaker can be used to amplify the receiver audio signal further or to relay the audio signal to a nearby room.

Advanced programming

By using advanced programming, you can customize the signal pattern from a specific transmitter and event, displaying the LED colour, sound and vibration pattern of your choice. The advanced programming overrides the radio key and pairs the units via the serial number. Please note that smoke alarms cannot be programmed for safety reasons.

Note: The transmitter must be activated as it is intended to be used in the system to generate the right signal. This means that you can't always use the transmitter test button (see **Default signal pattern** for the relevant transmitter).

Here is how you program the receiver:

- 1 Press and hold the test button on the receiver. The green and yellow Visit LEDs will start to blink alternately. While still holding down the button, activate the desired transmitter as intended. The power LED on the receiver will light up in yellow to show that you are in advanced programming mode. Release the button.
- 2 Scroll through the different Visit LED options by pressing the test button on the receiver. Select the desired Visit LED colour by holding down the test button until the power LED goes out and lights up again.
- **3** Scroll through the different **sound options** by pressing the test button on the receiver. Select the desired sound by holding down the test button until the power LED goes out and lights up again.
- 4 Scroll through the different vibration options by pressing the test button on the receiver (bed shaker required). Select the desired vibration pattern by holding down the test button until the power LED goes out and lights up again.
- 5 The receiver will now show the new Visit LED colour, sound and vibration pattern. Press the test button briefly to end the demonstration. After a short while, it will return to normal mode.

Deleting the advanced programming

Follow the procedure below to delete the advanced programming.

- 1 Hold down the test button on the receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- 2 Press the test button on the receiver 3 times in quick succession.
- 3 All Visit LEDs will light up for ~2 seconds to show that it has been deleted.

Troubleshooting

lf	Try this
The receiver seems to be turned off	 The batteries are depleted
The power LED is yellow when the receiver is activated.	 The battery level is low. Re
The receiver does not respond when a transmitter is activated, but works when I use the test button	 Check the transmitter batt Move the receiver closer t Check that the receiver is system, see Changing the
The receiver is activated for no apparent reason	 There is probably another Change the radio key on a
The receiver is too silent	 Turn up the volume using
The receiver is not flashing	 Check that the flash signal



ed. Replace them with 4 x 1.5V LR14 alkaline batteries.

Replace them with 4 x 1.5V LR14 alkaline batteries.

tteries and connections.

to the transmitter to make sure it's within radio range. s set to the same radio key as the other units in the Visit e radio key.


er Visit system installed nearby that triggers your system. all units, see Changing the radio key.

g the red volume dial on the top of the unit.

al switch on the back of the unit is set to the ON position.

Visit pager receiver

Buttons and controls



Belt clip Battery cover Battery eject Battery compartment + 0 \bigcirc \bigcirc BATTERY \bigcirc 6 1.5 V AAA ALKALINE BATTERY _

Technical specifications

In the box

- BE1470 Visit pager receiver
- Safety cord with clip
- Extra front label

Power and battery

- Mains power 7.5 V DC 200 mA via the charger
- Battery power 1.5 V AAA alkaline or 1.2 V AAA NiMH rechargeable battery
- Operation time Alkaline battery: 2 – 3 weeks NiMH battery: ~1 week
- Power consumption Active: ≤200 mA, Idle position: ≤1 mA

Optional safety cord with clip

Dimensions and weight

Height: 86 mm, 3.4"

Width: 57 mm, 2.2"

Depth: 29 mm, 1.1"

Visit LEDs

the following:

Weight: 70 g, 2.5 oz. incl. battery

The Visit LEDs normally indicates

Orange LED, pacifier symbol

Green LED, door symbol

Red LED, fire symbol

The baby monitor is activated

The door transmitter is activated

The phone transmitter is activated

Yellow LED, telephone symbol

The smoke alarm is activated

Environment

- For indoor use only Operating temperature 15° to 35° C, 59° to 95° F
- Relative humidity 5% to 95%, non-condensing

Frequency and coverage


- Radio frequency 314.91 MHz, 433.92 MHz or 868.30 MHz, depending on the region
- Coverage 50 - 250 m, 55 - 273 yd. depending on the radio frequency and the characteristics of the building


Accessories

- BE1260 Pager charger
- BE1270 Bed shaker

Getting started

- 1 Open the battery cover, fit the battery and close the cover again. Attach the pager to your belt using the belt clip. For extra security, use the supplied safety cord.
- 2 To test the radio link you need a Visit transmitter. Press the test button/s on the transmitter.
- 3 The pager starts to vibrate and lights up a Visit LED. If a bed shaker is connected during charging, it will vibrate. If nothing happens, see Troubleshooting.





Default signal pattern

When a transmitter is activated, the pager lights up a LED and starts to vibrate with a certain pace. This is called signal pattern. The transmitters determine the pattern, and the default is as follows:

Activated transmitter	Pager LED	Pager / bed shake
 Door transmitter 	Green	Slow ■□□□
Push button transmitter	Green	Slow ■□□□
 Telephone transmitter 	Yellow	Medium ∎□∎□
 Baby monitor 	Orange	Fast
Smoke alarm	Red	Long

Changing the signal pattern

The signal pattern can only be changed on the transmitters. See **Changing the signal pattern** for the relevant transmitter.

Changing the front label

If you want to use Visit for other purposes, the pager front label can be replaced with a customized one. Here is how it's done: • Open the battery cover, replace the original label with the supplied extra label and close the cover again.

Replacing the battery

When the battery LED starts to blink in yellow, the battery is nearly depleted. Here is how you replace it:

Open the battery cover and press the battery eject button to remove the old battery. Insert a 1.5 V AAA alkaline battery or a 1.2 V AAA NiMH rechargeable battery if you are using the BE1260 charger accessory.



er vibration

Visit pager receiver


Changing the radio key

If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference, you need to change the radio key on all units. The radio key switches are located on the *transmitters*.


Here is how you change the radio key:

- 1 Open the transmitter cover and move any radio key switch to the up = on position to change the radio key. See Changing the radio key for the relevant transmitter.
- 2 Press and hold the function button on the pager until the green and yellow Visit LEDs blink alternately. Release the button.
- 3 Press the test button/s on the transmitter within 30 seconds to send the new radio key.
- 4 All Visit LEDs on the pager blink 5 times to show that the radio key has been changed. It then returns to normal mode.

Please note: All Visit units must be set to the same radio key in order to operate as a group.







Of

Pager accessories


The pager can be complemented with the following accessories:

BE1260 Pager charger

Charges your pager during the night. Place it on the bedside table and connect up to two bed shakers. Please note that the pager will not vibrate when it's charging, but the Visit LEDs will act as usual.

BE1270 Bed shaker

Wakes you with vibrations if anything happens while you are asleep. Connect it to the pager charger and slide it under your pillow or mattress.





Warning! When using the pager charger ONLY USE RECHARGEABLE NIMH BATTERIES in the pager. Non-rechargeable batteries will start to leak if the pager is placed in the charger and the battery acid will damage the electronics. The resulting damage is not covered by warranty.

Advanced programming

By using advanced programming, you can customize the signal pattern from a specific transmitter and event, displaying the LED colour and vibration pattern of your choice. The advanced programming overrides the radio key and pairs the units via the serial number. Please note that smoke alarms cannot be programmed for safety reasons.

Note: The transmitter must be activated as it is intended to be used in the system to generate the right signal. This means that you can't always use the transmitter test button (see **Default signal pattern** for the relevant transmitter).

Here is how you program the pager:

- 1 Press and hold the function button on the pager. The green and yellow Visit LEDs will start to blink alternately. While still holding down the button, activate the desired transmitter as intended. The yellow battery LED on the pager will light up to indicate that you are in advanced programming mode. Release the button.
- 2 Scroll through the different Visit LED options by pressing the function button on the pager. Select the desired Visit LED pattern by holding down the function button until the battery LED goes out and lights up again.
- **3** Scroll through the different **vibration options** by pressing the function button on the pager. Select the desired vibration pattern by holding down the function button until the battery LED goes out and lights up again.
- 4 The pager will now show the new Visit LED colour and vibration pattern. Press the function button briefly to end the demonstration. After a short while, the pager will return to normal mode.

Deleting the advanced programming

Follow the procedure below to delete the advanced programming.

- 1 Hold down the function button on the pager until the green and yellow Visit LEDs blink alternately. Release the button. 2 Press the function button on the pager 3 times in quick succession.
- 3 All Visit LEDs will light up for ~2 seconds to show that it has been deleted.

Troubleshooting

Most problems with the pager can be solved quickly by following the advice below.


lf	Try this
The pager seems to be turned off	 The battery is depleted. F Important! If you NiMh battery in th
The battery LED blinks in yellow	 The battery level is low. R Important! If you NiMh battery in th
The pager does not respond when a transmitter is activated	 Check the batteries in the Move the pager closer to Check that the pager is se system, see Changing th
The pager is activated for no apparent reason	 There is probably anothe Change the radio key on



- Replace it with a 1.5V AAA alkaline battery. have a pager charger; only use a rechargeable 1.2 V AAA he pager.
- Replace it with a 1.5V AAA alkaline battery. have a pager charger; only use a rechargeable 1.2 V AAA he pager.
- ne transmitters.
- o the transmitter to make sure it's within radio range. set to the same radio key as the other units in the Visit he radio key.
- er Visit system installed nearby that triggers your system. all units, see Changing the radio key.

Visit wrist receiver

Buttons and controls



Technical specifications

In the box

- BE1560 Visit wrist receiver
- BE1570 Charger
- Elastic wrist band
- Power supply

Power and battery

- Mains power: 7.5 V DC/1500 mA
- Power consumption Receiver: Active: 100 mA, Idle: 3 mA Charger: Active: 650 mA, Idle: 70 mA
- Battery power Receiver: 1 x 1.2 V V40H rechargeable Charger: 4 x 1.2 V NiMH rechargeable
- Operating and charging time Receiver: ~30 h, Charging time: ~8 h Charger: Battery charging time: ~24 h

Dimensions and weight

- Receiver Charger Height 49 mm 100 mm
- Width 38 mm 95 mm 117 mm
- 12 mm Depth 185 g
- Weight 27 g

Visit LEDS

- The Visit LEDs normally indicates the following:
- Orange LED, pacifier symbol The baby monitor is activated
- Green LED, door symbol The door transmitter is activated
- Yellow LED, telephone symbol The phone transmitter is activated
- Red LED, fire symbol
- The smoke alarm is activated

Environment

- For indoor use only Operating temperature 15° to 35° C, 59° to 95° F
- Relative humidity 5% to 95%, non-condensing

Frequency and coverage


- Radio frequency 314.91 MHz, 433.92 MHz or 868.30 MHz, depending on the region
- Coverage 50 - 250 m, 55 - 273 yd. depending on the radio frequency and the characteristics of the building

Accessories

- BE1270 Bed shaker
- BE9086 External trigger cable

Getting started

- 1 Pull the battery tab on the charger and connect the power supply to the mains outlet. The mains LED lights up in green. Place the wrist receiver in the charger and charge it for at least 2 hours. The charging LED is green during charging.
- 2 To test the radio link you need a Visit transmitter. Press the test button/s on the transmitter.
- 3 The receiver starts to vibrate and lights up a Visit LED. If a bed shaker is connected during charging, it will vibrate. If nothing happens, see Troubleshooting.



Default signal pattern

When a transmitter is activated, the wrist receiver lights up a LED and starts to vibrate with a certain pace. This is called signal pattern. The transmitters determine the pattern, and the default is as follows:

Activated transmitter Wrist receiver LED Door transmitter Green

Push button transmitter Green Telephone transmitter Yellow Baby monitor Orange Smoke alarm Red

Changing the signal pattern

LED indications

When the wrist receiver battery is nearly depleted, the battery LED starts to blink in yellow. The charging time is up to 8 h. The charger is equipped with a battery backup and the charger LEDs indicate the following:

LED	Indication	Status
Charging LED	Green light	The receiver
Mains LED	Green light	The charger
Mains LED	Green blinks	The charger

39



Wrist receiver / bed shaker vibration

Slow	
Slow	
Medium	
Fast	
Long	

The signal pattern can only be changed on the transmitters. See Changing the signal pattern for the relevant transmitter.

- r battery is being charged.
- is powered by mains voltage.
- is powered by the battery backup.

Visit wrist receiver


Changing the radio key

If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference, you need to change the radio key on all units. The radio key switches are located on the *transmitters*.

Here is how you change the radio key:

- 1 Open the transmitter cover and move any radio key switch to the up = on position to change the radio key. See Changing the radio key for the relevant transmitter.
- 2 Press and hold the function button on the wrist receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- 3 Press the test button/s on the transmitter within 30 seconds to send the new radio key.
- 4 All Visit LEDs on the receiver blink 5 times to show that the radio key has been changed. It then returns to normal mode.

Please note: All Visit units must be set to the same radio key in order to operate as a group.




Accessories

The wrist receiver can be complemented with the following accessory:

BE1270 Bed shaker

Wakes you with vibrations if anything happens while you are asleep. Connect it to the charger and slide it under your pillow or mattress.

Note: The bed shaker only works when the wrist receiver is placed in the charger. The receiver will not vibrate during charging, but the Visit LEDs will act as usual.



Off

Press the bed shaker test button on top of the charger to try the bed shaker vibration.



Warning! The wrist receiver and charger can ONLY USE RECHARGEABLE NIMH BATTERIES. Nonrechargeable batteries will start to leak during charging and the battery acid will damage the electronics in the wrist receiver and charger. The resulting damage is not covered by warranty.

Advanced programming

By using advanced programming, you can customize the signal pattern from a specific transmitter and event, displaying the LED colour and vibration pattern of your choice. The advanced programming overrides the radio key and pairs the units via the serial number. Please note that smoke alarms cannot be programmed for safety reasons.

Note: The transmitter must be activated as it is intended to be used in the system to generate the right signal. This means that you can't always use the transmitter test button (see **Default signal pattern** for the relevant transmitter).

Here is how you program the wrist receiver:

- 1 Press and hold the function button on the receiver. The green and yellow Visit LEDs will start to blink alternately. While still holding down the button, activate the desired transmitter as intended. The battery LED on the receiver lights up in yellow to show that you are in advanced programming mode. Release the button.
- 2 Scroll through the different Visit LED options by pressing the function button on the receiver. Select the desired Visit LED colour by holding down the function button until the battery LED goes out and lights up again.
- 3 Scroll through the different vibration options by pressing the function button on the receiver. Select the desired vibration pattern by holding down the function button until the battery LED goes out and lights up again.
- 4 The wrist receiver will now show the new Visit LED colour and vibration pattern. Press the function button briefly to end the demonstration. After a short while, it will return to normal mode.

Deleting the advanced programming

Follow the procedure below to delete the advanced programming.

- 1 Hold down the function button on the receiver until the green and yellow Visit LEDs blink alternately. Release the button.
- 2 Press the function button on the receiver 3 times in guick succession.
- 3 All Visit LEDs will light up for ~2 seconds to show that it has been deleted.

Troubleshooting

Most problems with the wrist receiver can be solved quickly by following the advice below.

If	Try this
The receiver seems to be turned off	The battery is depleted.
The battery LED blinks in yellow	The battery level is low. C
The receiver is not charging	 Check that the receiver is is connected. The mains I Charge or replace the bac
The receiver does not respond when a transmitter is activated	 Check the transmitter bat Move the receiver closer Check that the receiver is system, see Changing the
The receiver is activated for no apparent reason	 There is probably anothe Change the radio key on
The bed shaker does not vibrate	 Check that the bed shake




- Charge or replace it with a VARTA V40H NiHM battery.
- Charge or replace it with a VARTA V40H NiHM battery.
- is placed correctly in the charger and that the power supply s LED and charging LED should be lit. ackup batteries with four 1.2 V NiMH rechargeable batteries.
- atteries and connections.
- r to the transmitter to make sure it's within radio range. is set to the same radio key as the other units in the Visit he radio key.
- er Visit system installed nearby that triggers your system. all units, see Changing the radio key.
- er is connected and that the receiver is placed in the charger.

Visit alarm clock receiver

Buttons and controls





Technical specifications

In the box

- BE1580 Visit alarm clock
- BE1272 Bed shaker with sound
- Power supply
- 4 x 1.2 V AAA NiMH batteries

Power and battery

- Mains power
 7.5 V DC / 1000 mA
 External power supply unit
- Backup batteries 4 x 1.2 V AAA NiMH rechargeable batteries
- Battery backup operating time
 ~ 24 h when fully charged
- Battery backup charging time ~ 10 h from fully depleted

Dimensions and weight

- Height: 108 mm, 4.3"
 Width: 121 mm, 4.7"
- width: 121 mm, 4.7
- Depth: 92 mm, 3.6"
- Weight: 390 g, 13.7 oz. incl. batteries

Visit LEDs

The Visit LEDs normally indicates the following:

- Orange LED, pacifier symbol The baby monitor is activated
- Green LED, door symbol The door transmitter is activated
- Yellow LED, telephone symbol The phone transmitter is activated
- Red LED, fire symbol The smoke alarm is activated

Output signals

3 - Landline telephone, 4 - Power supply

- Sound 100 dB @ 10 cm, 950 Hz – 3 kHz
- Four high-intensity flashing LEDs
- Bed shaker power: 2.0 4.0 VDC The bed shaker emits a sound

Frequency and coverage


- Radio frequency 314.91 MHz, 433.92 MHz or 868.30 MHz, depending on the region
- Coverage
 50 250 m, 55 273 yd. depending on the radio frequency and the characteristics of the building

Accessories

- BE1271 Bed shaker without sound
- BE9250 Mobile phone sensor

Getting started

- 1 Pull the battery tab and connect the power supply to the alarm clock and the mains outlet. Connect the bed shaker, tuck it under the pillow or mattress, and place the alarm clock on the bedside table.
- 2 To test the radio link you need a Visit transmitter. Press the test button/s on the transmitter.
- **3** The alarm clock lights up a Visit LED and starts to sound and flash. The bed shaker emits a sound and vibrates. If nothing happens, see **Troubleshooting**.



Default signal pattern

When a transmitter is activated, the alarm clock lights up a LED, sounds, flashes and the bed shaker starts to vibrate with a certain pace. This is called signal pattern. The transmitters determine the pattern, and the default is as follows:

Transmitter	Alarm clock			Bed shaker
Activated source	Visit LED	Sound	Flash	Vibration
Door transmitter / push button transmitter	Green	Door chime	Yes	Slow ■□□□
 Telephone transmitter / connected telephone 	Yellow	Ring signal	Yes	Medium ∎□∎□
 Baby monitor 	Orange	Baby melody	Yes	Fast IDIDID
Smoke alarm	Red	Fire horn	Yes	Long

Changing the signal pattern

The signal pattern can only be changed on the transmitters. See Changing the signal pattern for the relevant transmitter.

Settings

Flash and sound on/off

Press the flash and sound on/off button marked with \P/\P on the back of the alarm clock repeatedly to toggle between the options. A \mathscr{J} icon will appear on the clock face when the flash is turned off and a \mathscr{J} icon when the sound is muted.

Backlight intensity

Press the backlight button marked with 🚈 on the back of the alarm clock repeatedly to adjust the intensity in five steps.

Time format


Press the time format button marked with 12/24h on the back of the alarm clock to toggle between a 24h and a 12h setting.



Visit alarm clock receiver


Setting the time

Press the left dial and turn it to set hours. To set minutes, press and turn the dial again. Press once again to save your settings.



Activating the alarm

Press the alarm on/off button to activate the alarm. The alarm indicator lights up in red. To turn off the alarm, press the button again.



Alarm clock accessories

The alarm clock can be complemented with the following accessories:

BE9105 Telephone cord

Use it to connect the landline telephone to the alarm clock RJ11 input and be alerted when the telephone rings.

BE9250 Mobile phone sensor


Connect it to the ext. trig. input and place it on the display to be alerted by incoming calls or messages.

BE9024 Contact mat

Connect it to the ext. trig. input to be alerted when your spouse leaves the bed.


Setting the alarm


Press the **right** dial and turn it to set hours. To set minutes, press and turn the dial again. Press once again to save your settings.



Using the snooze and night light

Press the snooze button briefly to snooze the alarm (fire alarms cannot be snoozed for security reasons). Press and hold the snooze button for 3 seconds to turn on the night light. Press the button again to turn it off.






Changing the radio key

If your Visit system is activated for no reason, there is probably a nearby system that triggers yours. In order to avoid radio interference, you need to change the radio key on all units. The radio key switches are located on the transmitters.

Here is how you change the radio key:

- 1 Open the transmitter cover and move any radio key switch to the up = on position to change the radio key. See Changing the radio key for the relevant transmitter.
- 2 Press and hold the snooze button on the alarm clock until the green and yellow Visit LEDs blink alternately. Release the button.
- 3 Press the test button/s on the transmitter within 30 seconds to send the new radio key.
- 4 All Visit LEDs on the alarm clock blink 5 times to show that the radio key has been changed. It then returns to normal mode.

Please note: All Visit units must be set to the same radio key in order to operate as a group.




Troubleshooting

lf	Try this
The alarm clock seems to be turned off	 Check that the power si Charge the backup batt
The 🗇 symbol on the clock face starts to blink	 The power supply is dis Connect the power sup
A 🔏 symbol appears on the clock face	• The receiver detects no
The alarm clock does not respond when a transmitter is activated	 Check the transmitter b Move the alarm clock clip Check that the alarm clock that the alarm clock visit system, see Chang
The alarm clock is activated for no apparent reason	 There is probably anoth Change the radio key or
The alarm volume is too low	• The volume increases g








- supply is connected correctly. tteries for a couple of hours.
- sconnected and the backup batteries are nearly depleted. pply and charge the backup batteries for a couple of hours.
- backup batteries. Pull the battery tab, see Getting started.
- batteries and connections.
- closer to the transmitter to make sure it's within radio range. lock is set to the same radio key as the other units in the ging the radio key.
- her Visit system installed nearby that triggers your system. on all units, see Changing the radio key.
- gradually and reaches over 100 dB.

Pager charger

Buttons and controls




Technical specifications

In the box

- BE1260 Pager charger with pre-mounted backup batteries
- External power supply
- 1 x 1.2 V NiHM rechargeable battery intended for the pager
- Supporting legs and screws + plugs

Charging LED

- Green light: The pager is charging
- No light: The pager is fully charged or the power supply is not connected to mains power.



Dimensions and weight


- Height: 78 mm, 3"
- Width: 88 mm, 3.5"
- Depth: 43 mm, 1.7"
- Weight: 385 g, 16.3 oz. incl. batteries

Accessories

- BE1270 Bed shaker Connects up to two bed shakers
- BE9086 External trigger cable Use output B



- 2 Fit the supporting legs to the back of the charger and place it on a level surface. You can also mount it on the wall using
- 3 Connect the power supply to the mains outlet and place the pager in the charger. The charging LED is green during



Testing the connection

the supplied screws and plugs.

Note: Charge the backup batteries for 24 hours before using it with a bed shaker.

- 1 To test the radio link you need the pager and a Visit transmitter. Press the test button/s on the transmitter (see Testing the connection for the relevant transmitter).
- 2 If the pager is placed in the charger, it lights up a Visit LED and the bed shaker starts to vibrate. If nothing happens, see Troubleshooting.

Troubleshooting

Most problems with the charger can be solved quickly by following the advice below.

lf	Try this
The pager doesn't charge when it's placed in the charger	 Check that the pager will light up in green t
	 If the charging LED do Connect the power su
The bed shaker doesn't vibrate when the the pager is activated	 Check that the bed sh Check that the pager will light up in green to the pager should be address that the pager should be address to the
	 If the charging LED do Connect the power su



Warning! When using the pager charger ONLY USE RECHARGEABLE NIMH BATTERIES in the pager. Non-rechargeable batteries will start to leak if the pager is placed in the charger and the battery acid will damage the electronics. The resulting damage is not covered by warranty.

Power and battery

Battery power

Environment

For indoor use only

Mains power: 8 VDC / 800 mA

changed at a service center.

Normal charging time: ~ 6 h

Vibrator power: 2.0 – 4.0 VDC

Pager charging time

4 x 1.2 V NiMH rechargeable batteries

Note: The backup batteries must be

W. depleted backup batteries: ~ 24 h



1 Important! Replace the old alkaline battery with the supplied 1.2 V NiHM rechargeable battery.

charging and goes out when the pager is fully charged. Connect the bed shaker and tuck it under the pillow or mattress.

is positioned correctly in the charger. The charging LED to show that the pager is being charged.

loesn't light up, the backup batteries may be depleted. supply to mains power and charge the backup batteries.


haker is connected correctly to the charger.

is positioned correctly in the charger. The charging LED to show that the pager is being charged.

loesn't light up, the backup batteries may be depleted. supply to mains power and charge the backup batteries.

Mobile phone sensor

Buttons and connections



Model overview

BE9250 model

The BE9250 model has an adapter and connects to all Bellman & Symfon products with a 3.5 mm ext. trig. input.

Compatibility

- BE1580 Visit alarm clock
- BE1370 Pro alarm clock
- BE1350 Classic alarm clock

Technical specifications

- Battery power: 1 x AAA 1.5 V alkaline battery
- Optical detection: Activated when the display lights up
- Light sensitivity: Visible light >3 lux for longer than 2 s
- Connectors: Mobile phone sensor: 2.5 mm mono jack plug Adapter: 3.5 mm mono jack plug
- Cable length: 120 cm, 4'
- Sensor dimensions and weight: 24 x 34 x 24 mm, 20 g Adapter dimensions and weight: 53 x 25 x 18 mm, 27 g

BE9251 model

The BE9251 model connects to all Bellman & Symfon products that features a 2.5 mm ext. trig. input.

Compatibility


- BE1431 Telephone transmitter
- BE1441 Visit flash receiver
- BE1442 Visit flash receiver with battery backup
- BE1444 Visit flash receiver with battery backup and advanced smoke detector supervision

Technical specifications

- Placement: Horizontally on the display
- Optical detection: Activated when the display lights up
- Light sensitivity: Visible light >3 lux for longer than 2 s
- Connector: 2.5 mm mono jack plug
- Cable length: 120 cm, 4'
- Sensor dimensions and weight: 24 x 34 x 24 mm, 20 g


Using BE9250 with the alarm clock

- 1 Open the battery compartment, fit the supplied battery and connect the mobile phone sensor to the 2.5 mm input. 2 Connect the adapter to the 3.5 mm ext. trig. input on the back of the alarm clock.
- **3** Place the sensor on the mobile phone or tablet display and use e.g. a landline telephone to call the mobile phone. When the display lights up, the yellow Visit LED on the alarm clock blinks and it starts to sound, flash and vibrate.




Using BE9251 with the flash receiver

- 1 Connect the mobile phone sensor to the 2.5 mm ext. trig. input on the back of the receiver. Place it on the mobile phone or tablet display.
- 2 Use for instance the landline telephone to call the mobile phone. When the mobile phone display lights up, the yellow Visit LED on the receiver lights up and it starts to flash.




Using BE9251 with the telephone transmitter

- 1 Open the telephone transmitter front cover and connect the mobile phone sensor to the 2.5 mm ext. trig. input. Place it on the mobile phone or tablet display.
- 2 Use for instance the landline telephone to call the mobile phone. When the mobile phone display lights up, the transmitter top LED lights up in green to show that a radio signal is being transmitted.
- 3 The yellow Visit LED on the receiver lights up to show that the signal was received. In addition, it starts to sound, flash or vibrate depending on the receiver.










Bed shaker

BE1270

Connections



Technical specifications

Function

Wakes you with vibrations under the pillow or mattress. Requires no internal battery and connects to all Visit receivers and charger accessories.

In the box

BE1270 bed shaker

Power consumption

- Operating voltage 2.0 – 4.0 V DC from a Visit receiver
- Power consumption 250 – 750 mA

Cables and connectors

- Cable length: 2 m, 6.5'
- Connector: 3.5 mm mono jack plug

Dimensions and weight

- Height: 88 mm, 3.5"
- Width: 88 mm, 3.5"
- Depth: 27 mm, 1.1"
- Weight: 120 g, 4.2 oz.

Environmental requirements

- For indoor use only
- Temperature: 59° to 95° F, 15°-35° C
- Relative humidity: 5% -95% Non-condensing

Maintenance and care

- Clean with a dry cloth
- Do not use household cleaners, aerosol sprays, solvents, alcohol, ammonia, or abrasives

Compatibility

The bed shaker can be connected to the following Visit receivers and charger accessories:

Visit receivers

- BE1450 Visit portable receiver
- BE1441 Visit flash receiver
- BE1442 Visit flash receiver with battery backup
- BE1444 Visit flash receiver with battery backup and advanced smoke detector supervision
- BE1570 Visit wrist receiver charger

Accessories


 BE1260 Pager charger accessory for the BE1470 Visit pager receiver

Flash receiver wall bracket

Using a flash receiver

BE1441


- 1 Use the drilling template to mark and drill holes for the screws and plugs.
- 2 Remove the flash receiver table stand using the Allen key.
- 3 Fit the wall bracket on the wall. Attach the bottom and mount the flash receiver on the wall bracket.




Using a flash receiver with battery backup


BE1442/BE1444


- 1 Use the drilling template to mark and drill holes for the screws and plugs.
- 2 Remove the screw marked in red, located at the bottom of the flash receiver.
- 3 Fit the wall bracket on the wall. Re-attach the screw to mount the flash receiver on the wall bracket.











Visit accessories


Magnetic switch

Monitors the door and window

Mount the magnetic switch on the door or window frame and connect it to the telephone transmitter. When the magnets are separated, the transmitter signals the Visit receiver.

Technical specifications

- Dimensions 25 x 62 x 13 mm, 1" x 2.5" x 0.5"
- Weight 25 g, 0.9 oz.
- Connector 3.5 mm mono jack plug
- Cable length 0.5 m, 1.6'
- Contact breaker
 Open > 1 cm, 0.4" from the magnet
 Closed < 2 cm, 0.8" from the magnet
- Colour White
- Environment For indoor use only



Contact mat

Signals when someone steps on it

Place the contact mat by the front door or by the bed and connect it to the telephone transmitter or baby monitor to be alerted when someone enters a room or leaves the bed.

Technical specifications


- Dimensions 720 x 390 x 3 mm, 28" x 15" x 0.1"
- Weight 290 g, 10 oz.
- Connector 3.5 mm mono jack plug
- Cable length 200 cm, 6.6
- Dust proof and sealed to IP64 (not waterproof)

Operation

	Contact
--	---------

- Contact resistance 1 Ω (depending on pressure)
- Operating pressure Nominal 25 kg over 50 mm disc.
- Temperature range
- -10 to 70 °C, 14 to 158 °F

N/O normally open



Maximum ratings

- Contact rating 10 VASwitching voltage 25 VDC
- Switching current 0.25
- Carry current
- 0.25 Amps DC resistive 0.25 Amps DC resistive

External microphone

Extends the door transmitter reach

The external microphone can be used with the door transmitter when the sound source is located too far away from the internal microphone or when you need individual notifications from for example the doorbell and intercom.

Technical specifications

- Dimensions 33 x 36 x 4 mm, 1.3" x 1.4" x 0.2"
- Weight 15 g, 0.5 oz.
- Connector 3.5 mm mono jack plug
- Cable length BE9199: 2.5 m, 8.2
 - BE9200: 0.75 m, 2.5'
- Microphone type Piezoelectric
- Colour White
- Environment For indoor use only

External trig cable

Connects an external trigger source to Visit

The external trigger cable is used to connect an external trigger source to a Visit product. Use it for instance to connect an existing doorbell to the telephone transmitter and be alerted when someone rings the doorbell.


Technical specifications

- Weight 25 g, 0.9 oz.
- Connector 3.5 mm mono jack plug
- Cable length 0.5 m, 1.6
- Colour
 White

Voltage

- 🕖 2 30 VDC
 - 3 24 VAC

See the relevant Visit product section for detailed information.





For more than 25 years Bellman & Symfon of Sweden has been dedicated to improve the quality of life for people with hearing- and care related needs. Our people and partners are devoted to this mission and we work closely with healthcare professionals and leading experts to make an easy and independent living possible for everyone.

Bellman & Symfon Södra Långebergsgatan 30 421 32 Västra Frölunda Sweden Phone +46 31 68 28 20 Fax +46 31 68 28 90 Email info@bellman.com



 ${\rm \bigcirc}$ and TM 2016, Bellman & Symfon AB, all rights reserved. Bellman Visit and Bellman & Symfon are registered trademarks owned by Bellman & Symfon AB.