
# Bed shaker

BE1270

2	Čeština	7	Français	12	Norsk
3	Deutsch	8	Englsh	13	Polski
4	Dansk	9	Magyar	14	Português
5	Español	10	Italiano	15	Svenska
6	Suomi	11	Nederlands	16	Appendix

1

CE0402

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CZ

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
## User manual


GB

The bed shaker can be used with all the Visit receivers and chargers, as a means of waking up with help of vibrations.

### Getting started

- 1 Connect the bed shaker by inserting the connector in the bed shaker jack on the receiver or charger. (See below)
- 2 Press the test button on any of your Visit transmitters and check that the bed shaker vibrates.

**Note!** When you connect the bed shaker to the pager charger, the pager must be inserted in order for it to work.

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
## Technical information

APP

Operating voltage	2.0 - 4.0 V DC from a Visit receiver or charger
Power consumption	250 - 750 mA
Dimensions	88 x 88 x 27 mm (w x h x d)
Weight	120 g
Colour	White with red rubber anti-skid protector
Cable length	2 meters

Function	The BE1270 Bed shaker is a vibrator for indoor use that is placed in the bed to attract the attention of the user with vibrations.
Installation	The bed shaker is usually placed under a pillow or under the mattress. It is connected to a range of Bellman & Symfon products with vibrator outputs.
Testing	Connect the bed shaker to a Visit receiver or charger with vibrator output and follow the testing procedure.
Usage and storage	For indoor use only. Temperature during transport and storage: -10° to 50° C, 14° to 122° F. Temperature during operation: 0° to 35° C, 32° to 95° F. Relative humidity: 15% -90% (non condensing).
Care and cleaning	Maintenance free. Clean with a soft cloth. Do not use cleaners, solvents, ammonia or alcohol. Not reusable.

## Troubleshooting

If	Try this
The bed shaker does not vibrate	<ul style="list-style-type: none"><li>• Check that the plug is connected correctly.</li><li>• Check that the receiver is picking up signals.</li></ul>

## FCC compliance statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## Industry Canada statement

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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Revision: BE1270\_020MAN005  
Date of issue: 2020-10-28

TM and © 2020 Bellman & Symfon AB. All rights reserved.

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Bellman & Symfon Group AB