BE9026 Contact mat

- EN** Medical device information
- DE** Informationen zu Medizinprodukten
- FI** Lääketieteellisiä laitteita koskevat tiedot
- NL** Informatie over medische hulpmiddelen
- NO** Informasjon om medisinsk utstyr
- SE** Information om medicinteknisk utrustning

EN Medical device information

Read this first

Thank you for choosing a product from Bellman & Symfon – the world leader in alerting systems based in Gothenburg, Sweden. This leaflet contains important medical device information. Please read it carefully to make sure that you understand and get the best out of your Bellman & Symfon product. For more information about features and benefits, contact your hearing care professional.

About the BE9026 Contact mat

Intended purpose

This device is part of the Visit alerting system. The intended purpose of the system is to alert deaf and hard of hearing people of important signals in their home. The contact mat detects when someone steps on it. It typically activates a Visit transmitter that sends a signal to a Visit receiver which alerts the user with sound, flashes, or vibrations.

Intended user group

The intended user group consists of people of all ages with mild to severe hearing loss or deafness that need audio, visual or sensory amplification.

Intended user

The intended user is a person with mild to severe hearing loss or deafness.

Principle of operation





The contact mat should be placed on a flat, smooth surface. Connect it to the compatible Visit device (see Configuration options) to be alerted when someone steps on it.

 The Visit alerting system will not restore normal hearing and will not prevent or improve a hearing impairment or deafness resulting from organic conditions.

Important safety information

This section contains important information about safety, handling and operating conditions. Keep this leaflet for future use. If you are just installing the device, this leaflet must be given to the householder.

Hazard warnings

-  Failure to follow these safety instructions could result in fire, electric shock, or other injury or damage to the device or other property.
-  Keep this device out of reach of children under 3 years.
-  Do not use or store this device near any heat sources such as naked flames, radiators, ovens or other devices that produce heat.
-  Do not dismantle the device; there is a risk of electric shock.
-  Tampering with or dismantling the device will void warranty.
-  This device is designed for indoor use only. Do not expose the device to moisture.
-  Protect the device from shocks during storage and transport.
-  Do not make any changes or modifications to this device.
-  Use only original Bellman & Symfon accessories to avoid any electrical shock.
-  The contact mat should not be used in or relied on in safety related applications.
-  Cover the contact mat with a carpet or entrance mat to protect it from daily wear and to avoid trip or slip hazards.
-  Care should be taken not to stress the cable to pressure mat junction; this is especially important where long leads are fitted.

Information on product safety

-  Failure to follow these instructions could result in damage to the device and void the warranty.
-  Do not use the device in areas where electronic equipment is prohibited.
-  The device may only be repaired by an authorized service center.
-  If a serious incident occurs in relation to this device, contact the manufacturer and relevant authority.
-  If you encounter other problems with your device, contact the point of purchase, your local Bellman & Symfon office or the manufacturer. Visit bellman.com for contact info.

Operating conditions

Operate the device in a dry environment within the temperature and humidity limits stated in this leaflet. If the device gets wet or is exposed to moisture, it should no longer be regarded as reliable and should therefore be replaced.

Cleaning

Disconnect all cables before you clean your device. Use a soft, lint-free cloth. Avoid getting moisture in openings. Do not use household cleaners, aerosol sprays, solvents, alcohol, ammonia or abrasives. This device does not require sterilization.


Service and support

If the device appears to be damaged or doesn't function properly, follow the instructions in the user guide and this leaflet. If the product still doesn't function as intended, contact your local hearing care professional for information on service and warranty.

Warranty conditions

Bellman & Symfon guarantees this product for six (6) months from date of purchase against any defects that are due to faulty materials or workmanship. This guarantee only applies to normal conditions of use and service, and does not include damage resulting from accident, neglect, misuse, unauthorized dismantling, or contamination howsoever caused. This guarantee excludes incidental and consequential damage. Furthermore, the warranty does not cover Acts of God, such as fire, flood, hurricanes and tornadoes. This warranty gives you specific legal rights and you may also have other rights that vary with territory. Some countries or jurisdictions do not allow the limitation or exclusion of incidental or consequential damages, or limitations on how long an implied warranty lasts, so the above limitation may not apply to you. This guarantee is in addition to your statutory rights as a consumer. The above warranty may not be altered except in writing signed by both parties hereto.

Regulatory symbols

	With this symbol, Bellman & Symfon confirms that the product meets the Medical Device Regulation EU 2017/745.
	This symbol indicates the manufacturer's serial number so that a specific medical device can be identified. It's available on the product and gift box.
	This symbol indicates the manufacturer's catalogue number so that the medical device can be identified. It's available on the product and gift box.
	This symbol indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
	This symbol indicates that the user should consult this leaflet.
	This symbol indicates that it is important for the user to pay attention to the relevant warning notices in the user guides.
	This symbol indicates important information for handling and product safety.
	Temperature during transport and storage: -10° to 50° C, 14° to 122° F. Temperature during operation: -15° to 65° C, 5° to 149° F
	Humidity during transportation and storage: <90%, non-condensing. Humidity during operation: 15% – 90%, non-condensing
	Atmospheric pressure during operation, transportation and storage: 700hpa - 1060hpa
	This device is designed such that it functions without problems or restrictions if used as intended, unless otherwise noted in this leaflet.
	With this CE symbol, Bellman & Symfon confirms that the product meets EU standards for health, safety, and environmental protection as well as the Radio Equipment Directive 2014/53/EU.
	This symbol indicates that the product shall not be treated as household waste. Please hand over your old or unused product to the applicable collection point for the recycling of electrical and electronic equipment or bring your old product to your hearing care professional for appropriate disposal. By ensuring this product is disposed of correctly, you will help prevent potential negative effects on the environment and human health.

Technical specifications

- Weight:** 255g, 8.9 oz
- Dimensions:** 54 x 39 x 0.5 cm, 21" x 16" x 0.2"
- Standard operating pressure:** Nominal 5 kg over 50 mm disc.
- Colour:** Black
- Cable length:** 2 meters, 6.6 ft
- Connector:** 3.5 mm mono plug
- In the box**
 - BE9026 Contact mat

Configuration options

This contact mat can be configured with the following transmitters and accessories:

- | | |
|--|--|
| Compatible transmitters <ul style="list-style-type: none"> ■ BE1431 Visit telephone transmitter ■ BE1432 Visit mobile phone transmitter ■ BE1433 Visit mobile phone transceiver ■ BE1491 Visit baby monitor ■ BE1370 Alarm clock pro ■ BE1580 Visit Alarm clock | Compatible accessories <ul style="list-style-type: none"> ■ BE1510 Repeater ■ BE1520 Bridge |
|--|--|

 For detailed product information, see the corresponding user manual.

Compliance information

Hereby Bellman & Symfon declares that, in Europe, this product is in compliance with the essential requirements of the Medical Device Regulation EU 2017/745 as well as the directives and regulations listed below. The full text of the declaration of conformity can be obtained from Bellman & Symfon or your local Bellman & Symfon representative. Visit bellman.com for contact information.

Medical Device Regulation (MDR)
Electromagnetic Compatibility Directive (EMC)
Restriction of Hazardous Substances Directive (RoHS)
REACH Regulation
Waste Electrical & Electronic Equipment (WEEE)

ISO Certification of legal manufacturer

Bellman is certified in accordance with SS-EN ISO 9001 and SS-EN ISO 13485.
SS-EN ISO 9001 Certification Number: CN19/42071
SS-EN ISO 13485 Certification Number: CN19/42070

Certification Body

SGS United Kingdom Ltd
Rossmore Business Park Ellesmere Port Cheshire CH65 3EN UK

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All rights reserved.
Bellman & Symfon Group AB
© 2021 Bellman & Symfon AB.