BELLMAN AUDIO DOMINO

EN

DE

DK

ES

FR

JP

PL

SE

SI-CN

TR-CN

APP

PERSONAL HEARING SYSTEM

## Innehåll

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
# Appendix

## Settings

Bellman Audio Domino has more settings than those you can access in user mode. By going into the Settings Menu you can adjust settings that you do not normally need to change or which are adjusted when first used and never need to be changed again.

### Settings menu

Start the Receiver Unit with the On/Off button  $\odot$  (6) for about 3 seconds.

Hold down the On/Off button  $\odot$  (6) and press the Volume Button — (10) within 1 second for at least 3 seconds.


You are now in Balance Control mode.

### Overview of the Settings Menu

In the settings menu you can move between the different settings by

pressing the **RF** button (7) or the **MIC** button (5).

Each press on the **RF** button moves you clockwise on to the next setting in the list and the **MIC** button takes you anti-clockwise to the previous setting in the list.


When you have finished adjusting your settings, give a short press on the On/Off button  $\odot$  (6) to return to normal user mode.

## Balance control

The balance between the right and left-hand channels can be set to compensate, for example, for a hearing loss in one ear.


Press the **<** button (11) to adjust the balance on the left and the **>** button (12) to adjust the balance on the right.

The standard setting is centred balance.

## Basic attenuation

Different headphones have different sensitivity levels and you can compensate for this with basic attenuation. What you actually set is the attenuation you require at the maximum power output you want to have (which is therefore 0 dB attenuation). This settings is useful when you want to adjust the maximum volume in user mode for a pair of headphones or in-ear phones with a higher sensitivity than the BE9122 Bellman Stereo Headphones or, for example, if you want to limit the volume for a child.


MAX **<**

Use the **>** button (12) to reduce the basic attenuation (increase the volume) and the **<** button (11) to increase the basic attenuation (reduce the volume). The basic attenuation is adjusted in increments of 3 dB.

The standard setting is 0 dB attenuation, which corresponds to the maximum output signal.

« STEREO »


## Stereo/Mono

The line-in signal can be rendered in either Stereo or in Mono. The speech intelligibility can be increased in certain listening situations such as when watching TV if the signal is in Mono.

Press the **+** button to select Stereo.

MONO **>**


Press the **—** button to select Mono.

The standard setting is Stereo.

## Bass reduction


Low-frequency interference such as fans and road noise can be tiring if you are exposed to it for a long period. By activating bass reduction, low-frequency interference and sounds are reduced.

BASS 

Press the **+** button (9) to activate Bass Reduction.

BASS 


Press the **-** button (10) to deactivate Bass Reduction and get level frequency response.

The standard setting is deactivated Bass Reduction.

## Adjusting the Tone Control


The Tone Control Blocking function has been developed in case you want to make it easier for the user to use Bellman Audio Domino. This setting can be used to choose whether the user should be able to adjust the tone control or not. If you choose to switch off the tone control buttons, you can preset the tone control that will be used in user mode.

### Activating and deactivating Tone Control in User Mode

Press the **+** button (9) to activate the tone control buttons and the **-** button (10) to deactivate them.

The standard setting is tone control activated.


### Presetting the Tone Control

Press the **<** button (11) and the **>** button (12) to preset the tone control to the level that will be used in user mode.

The standard setting is Tone Control centred (level frequency response).


Please note that if the tone control is on, the preset value will be changed as soon as the Tone Control is changed in user mode.

### Resetting to the factory setting

You can use this to reset all settings to the standard settings that the products were supplied with from the factory.

Press and hold down the **+** button (9) until all the LEDs on the display are lit. After resetting is completed, the display will blink for a short period.



### Version information

The product version information is presented here via a pattern of lit and blinking LEDs. This internal information is only intended as an aid for servicing and support.

### Pairing of units

The Transmitter and Receiver Units in the Bellman Audio Domino system must have the radio link established in order for radio communication to operate. To make it easier to get started, a Receiver Unit and a Transmitter Unit which are paired together are always supplied so that they will work together. Where you want to pair a Receiver Unit with a Transmitter Unit which is different from the Receiver Unit it was originally supplied with, you can do this in the Settings Menu. Proceed as follows to pair up the units:

- Start both the units with the On/Off buttons **⏻** (6, 18).
- Go into the Settings Menu.
- Hold down the **CONNECT** button (19) on the Transmitter Unit.
- Hold down both the **MIC** button (5) and the **RF** button (7) on the Receiver Unit.

A "running light" on the Receiver Unit's Volume and Tone LED Indicator (8) moving towards the centre indicates that the link has been established.

APP


Please note! If you pair up a Transmitter Unit and a Receiver Unit in the Settings Menu, this replaces the pairing carried out in the factory. This can be useful when changing a Transmitter Unit under guarantee or if you want to pair up several Receiver Units to a Transmitter Unit (e.g. in a conference room).

## Educational Mode

The Domino Pro has a mode especially developed for use in an educational setting. This mode sets the microphone in the transmitter to be active by default.

The microphone in the receiver can be momentarily activated to allow conversation with peers and support staff around the user.

- By pressing and holding the MIC button, the sound from the local microphone will be mixed with the sound from the transmitter microphone
- When the MIC button is released, the receiver will go back to the transmitter microphone again

### Turn On the Educational mode

- 1) Enter the “Settings menu” (refer to Page 276)
- 2) Use the MIC (5) button to move between settings until the LED's over the MIC (5) and RF (7) buttons flash (this is the “Resetting to the factory setting” menu).
- 3) Press and hold the Volume down — button (10) for 20 seconds.
- 4) The LED lightbar will flash to acknowledge the change of mode.
- 5) Press the power button once shortly to set education mode and return to the normal menu.

### Turn Off the Educational mode

The Educational mode is turned OFF by restoring Factory settings (refer to Page 282).

Note: The factory default setting is Educational mode OFF.

## Testing

It is easy to test the Bellman Audio Domino Personal Hearing System. If Bellman Audio Domino does not operate as described below, you can carry out further troubleshooting as instructed under Troubleshooting.

### How to test

The testing of Bellman Audio Domino has been divided into three parts. One part for the Receiver Unit for which only the Receiver Unit is required, one part for the Transmitter Unit for which both units are required and finally one for charging the two units.

### Receiver Unit

- Switch on the unit and make sure that the Battery Status LED (3) lights up green.
- Connect the headphones to the Headphone Socket (13) and scrape carefully on the round microphone (1) with your finger and check that there is a scraping sound in the headphones.
- Press the **MIC** button (5) and make sure that the Microphone LED (2) changes colour and lights up orange.
- Scrape carefully on the small microphone hole on the top (1) and the round microphone (1) with your finger and check that there is a scraping sound in the headphones.
- Press the buttons (9, 10, 11, 12) and check that the Volume and Tone Indicator changes.

## Transmitter Unit

- Switch on the unit and make sure that the Battery Status LED (17) lights up green.
- Press the **RF** button (7) and make sure that the External Sound LED (4) lights up green.
- Carefully scrape on the round microphone hole (16) on the top of the Transmitter Unit with the finger and check that there is a scraping sound in the headphones.
- Press the **RF** button (7) and make sure that the External Sound LED (4) changes colour and lights up orange.
- Carefully scrape on the round microphone (16) on the Transmitter Unit with the finger and check that there is a scraping sound in the headphones.
- Connect an external sound source to the External Sound Source Input (20) and check that the sound from the sound source is heard in the headphones. The external sound source input (20) can be tested, for example, by connecting a CD player or a radio to the input via the audio cable supplied. The transmitter unit automatically switches to an external sound source when a connector is connected to the External Sound Source Input (20).

## Charging

- Connect the charger to the mains socket and check that the LED on the charger light up green.
- Connect the charging cables to both the charger and to the charging sockets on the Receiver Unit and the Transmitter Unit (14, 21), see the fold-out section of the manual.
- Check that the LED on the charger change from green to red.
- Check that the Battery Status LEDs on both the Transmitter Unit and Receiver Unit light up orange.

## Technical information

### Power supply

Battery power	Built-in rechargeable lithium-ion battery 3.7 V, 860 mAh
	 <b>CAUTION:</b> <i>Danger of explosion if battery is incorrectly replaced. Replace only with the same or equivalent type.</i>
Operating time, fully charged	<ul style="list-style-type: none"><li>• Up to 12 hours using BE9124 earphones</li><li>• Up to 8 hours using BE9159 neckloop</li></ul>
Transmitter charging time	4 hours
Receiver charging time	4 hours

*Please note! The battery capacity may deteriorate over time. If you find that the operating time is too short, please contact your service center for replacement of the rechargeable battery.*

### Input signals

Microphone	Internal electret condenser type (omni-directional or directional)
External sound source input	2.5 mm stereo tele jack plug 0 dBV sensitivity (max input level)  Automatically selects Aux. in or Ext. mic.

### Output signals

Headphone socket	3.5 mm stereo jack plug
------------------	-------------------------

Output power	100 mW @ 16 Ω load
Distortion	BE9122 Headphones 0.6 % THD BE9124 Ear Phones 0.08 % THD
Dynamic range	95 dB SNR (electrical)
Frequency range	40 – 10.000 Hz (electrical)
Max output level, saturation	
BE9122 Headphones	BE2210: 124.3 dB SPL / BE2230: 124.3 dB SPL (@ 2.5 kHz)
BE9124 Ear Phones	BE2210: 132.0 dB SPL / BE2230: 130.5 dB SPL (@ 5 kHz)
Full on gain	(60 dB SPL input signal @ 1 kHz):
BE9122 Headphones	BE2210: 42.5 dB / BE2230: 44.0 dB
BE9124 Ear Phones	BE2210: 55.0 dB / BE2230: 57.5 dB
Full on gain	(60 dB SPL insignal, HFA):
BE9122 Headphones	BE2210: 46.5 dB / BE2230: 49.0 dB
BE9124 Ear Phones	BE2210: 51.5 dB / BE2230: 56.0 dB

### Radio function

Digital radio technology	2.4 GHz frequency band, three channels Co-exists with other equipment Automatic channel selection 128 bits encryption
Topology	Point-to-point or Point-to-multipoint
Diversity	2 antennas for reception, one antenna for transmission (Rx + Rx/Tx)
Transmitting output power	Max. 25 mW

## Radio function

Operating range radio	Up to 25 meters indoors clear line of sight, up to 100 meters outdoors clear line of sight. Walls, concrete and metal reduce the range.
The following may also affect the performance:	More than one Receiver is connected to a Transmitter and at least one of the Receivers is too far away from the Transmitter. Too many Receivers are connected to the same Transmitter. Too many transmitters on 2.4 GHz or too many Domino systems close to one another.

## Features

Volume control	9 steps (5 dB/step)
Tone control	+/- 10 dB in 9 steps, cutoff frequency 1,6 kHz
Power switch	On/Off toggle button
Microphone and RF-mode	Separate mode select buttons
R/L-channel balance setting	0-20 dB programmable attenuation (in steps of 5 dB)
Maximum volume setting	0-21 dB programmable attenuation (in steps of 3 dB)
Bass reduction setting	Attenuation of frequencies below 400 Hz with 10 dB

Dynamic compressor:	10 channels Dynamic gain range 0-35 dB Compression 3:1 Expansion 1:1,5 Attack time 20 ms <3kHz, 25 ms ≥3kHz Release time 110 ms <3kHz, 200 ms ≥3kHz
Dynamic noise reduction	10-channel adaptive noise reduction
Feedback cancellation	Adaptive phase shift type
Digital signal processing	22.050 kHz sampling frequency, 16-bit resolution in stereo
Additional information	Stereo/Mono setting Reset of factory default settings High-polished gloss coating for reduction of handling noise



## Additional information

For indoor use only.

Dimensions H x W x D	99 x 48 x 22 mm
Weight	BE2210: 80g, BE2230: 80g
Colour	Black/dark grey metallic, silver-coloured buttons.

## Measurements

Frequency response	Acoustical (BE9122 Headphones, and BE9124 Ear Phones. See diagrams on next page.
--------------------	---

## Troubleshooting

Problem	Solution
The sound is louder in one ear than in the other.	Adjust the Balance in the Settings Menu.
The sound volume is low even though the volume is set to max.	Adjust the Basic Attenuation in the Settings Menu.
The sound is only in Mono even though it is connected via a stereo cable to a sound source which provides stereo sound.	Change to Stereo in the Settings Menu.
The bass reproduction is weak.	Switch off Bass Reduction in the Settings Menu.
The Tone Control cannot be adjusted.	Use the Settings Menu to activate the Tone Control buttons or to adjust the sound.
The Receiver Unit is behaving strangely.	Reset the unit to the factory settings in the Settings Menu.
The Receiver Unit loses its link to the Transmitter Unit when you switch off the Receiver Unit.	Pair up the Receiver Unit and Transmitter unit permanently in the Settings Menu.
Humming occurs when connecting the stereo cable during charging.	Only use the charger supplied, BE9142.

## FCC compliance statement

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- *Reorient or relocate the receiving antenna.*
- *Increase the separation between the equipment and receiver.*
- *Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.*
- *Consult the dealer or an experienced radio/TV technician for help*

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

### Statement for Industry Canada

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

(1) This device may not cause interference.

(2) This device must accept any interference, including interference that may cause undesired operation of the device.

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This equipment complies with IC exposure limits set forth for an uncontrolled environment. The device can be used in portable exposure conditions. This equipment shall be installed and operated with minimum distance 20cm between the radiator & body.

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