

# Stetoclips

BE9125



## Features

**High amplification**  
Up to 140 dB sound level

**User friendly**  
Easy to put on and take off

**Lightweight**  
Weigh only 38 g (1.3 oz)

## For extreme sound amplification.

The stetoclips are specially selected for the Bellman & Symfon range of digital speech amplifiers and listening systems. Designed especially for hard of hearing people, they are capable of producing a sound level in mono of up to 140dB.

The stetoclips are lightweight, comfortable to wear and very easy to put on and take off.

### Technical specifications

■ Frequency range	800 – 3000 Hz
■ Output power	124 dB SPL @ 1kHz, 50 mW
■ Connector	3.5 mm stereo jack plug
■ Cable length	1.5 m, 4.9'
■ Weight	38 g, 1.3 oz
■ Colour	Black

### Caution

The stetoclips are designed especially for hard of hearing people. They can produce high sound levels which could adversely affect your hearing if used inappropriately.

© 2016, Bellman & Symfon AB. All rights reserved.  
Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB. This leaflet is subject to change without prior notice.

BE9125\_001DAS002

[bellman.com](http://bellman.com)