

Getting Started

What you'll need

A smartphone with iOS 10 /
Android 4.4 or later.



What's in the box

- | | |
|-------------------------|------------------|
| 1 Telephone transmitter | 3 Screw and plug |
| 2 Telephone adapter | 4 Telephone cord |

1



2



3



4




* The appearance of your telephone adapter may differ.
© and TM 2018 Bellman & Symfon AB. All rights reserved.

Ready. Steady. Go.

1 Turn on Bluetooth®


On your iOS device

Swipe up from the bottom of the screen and tap the Bluetooth icon.



On your Android device

Swipe down from the top of the screen and tap the Bluetooth icon.



2 Set up your transmitter

Download the Visit app from the App Store® or Google Play™ and follow the simple setup instructions.



Download on the
App Store



GET IT ON
Google Play