

Vibio Bed Shaker

BE1221-P02



Features

Powerful vibrations

Wakes even heavy sleepers

Bluetooth®

Easy-to-use app included

Smart notifications

Alerts for calls and messages

Snooze strap

Snooze or stop with a gentle pull

Set and forget

Works even without the phone

A great way to start your day.

Whether you are ready to drop after the daily grind or just need a catnap to recharge, Vibio will make sure you get up on time. Set the alarm in the easy-to-use app, available for both iOS and Android, and tuck Vibio under your pillow or mattress. When it's time to start your day, Vibio will wake you with silent but powerful vibrations. Vibio is the perfect choice for heavy sleepers, couples and people with hearing loss.

Vibio can even notify for calls and messages during the nap, so you don't miss out on anything important. Not ready to get out of bed yet? Pull the snooze strap under the pillow for a few more minutes in the sack.



MD

CE-marked medical device under EU MDR 2017/745

Technical specifications

Dimensions and weight

Bed Shaker

Ø: 94 × 94 mm, 3.7 × 3.7"

Height: 27 mm, 1.1"


Weight: 153 g, 5.4 oz.



Snooze strap

Length: 120 mm, 4.7"

Buttons and indicators



Compatibility

- iOS15 / Android 8 and later
- iPhone 6S and later
- Modern Android phones

Frequency and coverage

- Bluetooth® frequency: 2400 – 2483.5 MHz
- Bluetooth® coverage: Up to 30 m, 33 yd, free field

Environmental requirements

- For indoor use only
- Temperature: 0° to 35° C, 32° to 95° F
- Relative humidity: 15% to 90% (non-condensing)

Cables and connections

- Cable length: 1 m, 39"
- Contacts: USB type A to USB Micro-B

Power and battery

- Battery power: 1 × 3.7 V AA Lithium-ion
- Power consumption: < 1000 mA
Idle: < 100 µA
- Operating time: Up to 40 days / charge
- Charging time: Up to 1.5 hour @ 500mA
- Charging current: Up to 500 mA, 5V DC

Maintenance and cleaning

- Maintenance free. Clean with a dry cloth.
- Do not use window cleaners, household cleaners, aerosol sprays, solvents, alcohol, ammonia, or abrasives.

Connectivity

- Bluetooth® 5

Approvals



This product is a medical device in conformity with Regulation (EU) 2017/745 (Medical Device Regulation).

In the box

- BE1221-P02 Vibio Bed Shaker
- USB charging cable



© 2025 and ™ Bellman & Symfon AB

Bellman and Bellman & Symfon are registered trademarks owned by Bellman & Symfon AB, all rights reserved. This leaflet is subject to change without prior notice.

BE1221-P02_033DSH1.0_EN

bellman.com