


Visit Pager Charger

BE1260-C



Features

Night-time alerts
Alerts you with vibrations

Bed shaker outlets
Supports two bed shakers

Battery backup
Works even during power loss



Charger for the Visit Pager.

The Visit Pager Charger from Bellman & Symfon is specially designed for the Visit System and is used to recharge the Visit Pager battery.

The Pager alerts you even during charging, just connect a bed shaker to the charger to wake up to vibrations if a Visit transmitter is activated while you are resting. The charger features backup batteries and remains fully functional during power loss, so you can get a peaceful sleep.

Technical specifications

- Dimensions 78 x 88 x 43 mm, 3.1" x 3.5" x 1.7"
- Weight, incl batteries and power adapter 200 g, 7 oz.
- Mains power 7.5 VDC / 1500 mA
- Backup battery Internal NiMH batteries
- Vibrator power 2.0 - 4.0 VDC
- Accessories:
 - BE1470-868 Visit Pager Receiver
 - BE1270-P02 Bed Shaker Connect up to two
- For indoor use only

Approvals

MD This product is a medical device in conformity with Regulation (EU) 2017/745 (Medical Device Regulation).



© 2025 , Bellman & Symfon AB. All rights reserved.
Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB. This leaflet is subject to change without prior notice.

BE1260-C_001DSH1.0_EN bellman.com