

Watch Receiver

BE3330



Your home on your wrist.

This Watch Receiver combines discreet vibrations with clear notifications in a sleek, lightweight design. It lasts up to one week on a single charge and offers an impressive range of up to 250 meters in open field conditions.

Whether you are in the basement or enjoying a cozy dinner in the garden, you will stay connected and within reach.

The Watch Receiver is made from lightweight materials that withstand everyday wear and tear, are gentle on the skin and comfortable to wear throughout the day. It is packed with useful features such as an alarm clock, timer and call-for-attention function. Simply hold down the back key and the Visit receivers around your home will sound, flash or vibrate to notify your loved ones.

Bellman.



Features

Discreet vibrations
Silent and non-intrusive

Clear signal identification
Symbols and vibration patterns

Comfortable to wear
Materials that are kind to the skin

Easy to customize
Exchangeable bands and dials

Operating time
Up to a week on a single charge

MD

CE-marked medical device
under EU MDR 2017/745



Technical specifications



Dimensions and weight

Lug width	50 mm
Width	45 mm
Depth	10 mm
Weight	43 g (Watch), 14 g (Charger)
Watch material	ABS + PC
Band material	Silicone



Overview



Features

- **Compatible with Visit units**
Via the BE1521-868/-433/-315 Bridge transceiver
- **Clear alert identification**
Dynamic symbols and vibration patterns
- **Multiple clock faces**
Both analogue and digital variants
- **Out of range alert**
Alerts you if connection is lost
- **Do Not Disturb**
Useful, when you want to take a nap
- **Night Mode**
Reduces the eye strain during the dark hours
- **Call for attention**
Press a button to call for your family's attention
- **Safe Mode**
When functional minimalism is key
- **Splash proof (IP67)**

Power and battery

- Lithium polymer rechargeable battery
3.7 V, 200 mAh
- Operating time
Up to one week
- Charging time
Up to two hours

Communication and coverage

- Bluetooth frequency
2400 – 2483.5 MHz
- Bluetooth transmission power
3mW (max.)
- Bluetooth range
Up to 250 m (270 yd)

In the box

- BE3330 Watch Receiver
- BE9292 Pocket Charger
- BE9419 Wrist Band



Approvals

MD This product is a medical device in conformity with Regulation (EU) 2017/745 (Medical Device Regulation).

Environmental requirements

- Operating temperature: 0° to 35° C, (32° to 95° F)
- Relative humidity: 15% to 90%
(non condensing)

© and TM 2026 Bellman & Symfon AB
Bellman and Bellman & Symfon are registered trademarks owned by Bellman & Symfon AB, all rights reserved. This leaflet is subject to change without prior notice.