

Microset™

BE9129



Features

Hands-free listening
Tuck away Mino in your pocket

Small and lightweight
For a discreet look

Works out of the box
No fitting procedure required

All-in-one microphone and earbud for the Mino.

The Microset is specially designed for the Mino speech amplifier. Put the Mino in your inner pocket or by using a belt clip, and use the lapel microphone to pick up the sound. It includes an omnidirectional lapel microphone and a single earphone.

Earphone specifications

- Sensitivity 113(± 3) dB/mW (@ 1kHz)
- Frequency response 50 Hz to 20 kHz
- Maximum power 10 mW
- Impedance 32 ohm

Microphone specifications

- Type Electret condenser microphone
- Directivity Omnidirectional
- Sensitivity -47 (±3) dB1
- Frequency response 110 Hz to 16 kHz
- S/N ratio >58 dB
- Bias voltage Max. 10 VDC

Cable and connector

- Cable length 95 cm, 3.1'
- Connector 3.5 mm stereo plug

Compliance

This accessory is designed for use with Bellman & Symfon systems, which are compliant with applicable FCC and ISCED Canada regulations.



© 2025, Bellman & Symfon AB. All rights reserved.
Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB. This leaflet is subject to change without prior notice.

BE9129_001DSH1.0-EN

bellman.com