External microphone

BE9136







Features

Use it during meetings To improve speech intelligibility

Use it for TV-listening

For clear and distinct sound

Adjust the cable length

With optional extension cable

For meetings and TV-listening.

The external microphone is specially developed for the Maxi and Mino digital speech amplifiers. It can be used as a remote microphone during meetings or conferences to improve the speech intelligibility in busy environments.

The external microphone can also be used for listening to the TV. Just attach the microphone close to the TV's loudspeaker and enjoy the loud and clear TV-sound.

Technical specifications

Cable length 1 m, 3.2'
 Connectors 2.5 mm mono jack plug

Directivity Omnidirection

Sensitivity -44dB

■ S/N >58dB


Weight 12 g, 0.4 oz.

Colour Black

■ In the package

Tie-clip Microphone foam cap

External microphone



Optional extension cable

BE9408 Cable for tie-clip microphone

Cable length 4 m, 13.1

Connectors2.5 mmmono jack plug

© 2020, Be Ilman & Symfon AB. All rights reserved.

Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB. This leaflet is subject to change without prior notice.

BE9136_003DSH001

bellman.com