

# External Microphone

BE9136



### Features

**Use it during meetings**  
To improve speech intelligibility

**Use it for TV-listening**  
For clear and distinct sound

**Adjust the cable length**  
With optional extension cable


## For meetings and TV-listening.

The External Microphone is specially developed for the Maxi, Mino and Domino digital speech amplifiers. It can be used as a remote microphone during meetings or conferences to improve the speech intelligibility in busy environments.

The External Microphone can also be used for listening to the TV. Just attach the microphone close to the TV's loudspeaker and enjoy the loud and clear TV-sound. Extension cable available as accessory for getting closer to the source.

### Technical specifications

- Cable length 1 m, 3.2'
- Connectors 2.5 mm mono plug
- Directivity Omnidirectional
- Sensitivity -44dB
- S/N >58dB
- Weight 12 g, 0.4 oz.
- Colour Black
- Accessories BE9408 Extension cable for Microphone.  
Cable length 4 m, 13.1'
- In the package BE1936-P02 External Microphone  
Tie-clip  
Microphone foam cap



### Compliance

This accessory is designed for use with Bellman & Symfon systems, which are compliant with applicable FCC and ISED Canada regulations.

© 2025, Bellman & Symfon AB. All rights reserved.  
Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB. This leaflet is subject to change without prior notice.