

In-Ear Headset

BE9229-P02



Features

Remote-button
With an integrated microphone

Closed drivers
Minimizes the risk of feedback

Small and lightweight
For a discrete look



Talk and listen undisturbed.

This In-ear Headset is specially designed for Maxi Pro and can be used to listen to music or talk on the phone undisturbed. It features a remote-button with a fully integrated microphone and offers high-clarity stereo sound with great detail.

The headset is lightweight and provide excellent isolation from ambient noise, increasing the overall speech intelligibility.

Technical specifications

- Frequency range20 Hz - 20 kHz
- Impedance16 Ω +/- 10%
- Sensitivity93 dB/mW @ 1 kHz SPL
- Connector3.5 mm stereo plug
- Cable length1.25 m, 49"
- Weight16 g, 0.6 oz
- ColourBlack
- Tangle-free cable
- Two sizes of ear tips included

Approvals

MD This product is a medical device in conformity with Regulation (EU) 2017/745 (Medical Device Regulation).

© 2025, Bellman & Symfon AB. All rights reserved.
Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB. This leaflet is subject to change without prior notice.

BE9229-P02_001DSH1.0 bellman.com