# **Stetoclips**









**Features** 

**High amplification** Up to 145 dB sound level

**User friendly** Easy to put on and take off

Lightweight Weigh only 40 g (1.4 oz)

## For extreme sound amplification.

The stetoclips are specially selected for the Bellman & Symfon range of digital speech amplifiers and listening systems. Designed especially for hard of hearing people, they are capable of producing a sound level in mono of up to 145dB.

The stetoclips are lightweight, comfortable to wear and very easy to put on and take off.

#### **Technical specifications**

- Frequency range
- Output power 145 dB SPL @ 1kHz, 50 mW
- Connector 3.5 mm stereo jack plug
- Cable length
- Weight
- Colour

### Caution

20 – 5000 Hz

0.8 m, 2.7 ft

40 g, 1.4 oz

Black

The stetoclips are designed especially for hard of hearing people. They can produce high sound levels which could adversely affect your hearing if used inappropriately.

#### © 2020, Bellman & Symfon AB. All rights reserved.

Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB. This leaflet is subject to change without prior notice.

BE9403\_002DSH001