


The alarm clock from Bellman & Symfon

The strong, silent type



Extra loud, audible alarm

Ascending multiple frequency signal that delivers up to 100+dB



High intensity flash

Four bright flashing LED lights that are effective even during the day



Powerful bed shaker

Vibrates and emits a sound that can be heard through the pillow



Night light beacon

Mild blue light that guides you back to bed if you are up at night



Safety backup power

Keeps all the functions, even the bedshaker, working during power loss



Smart snooze

Automatically decreases the snooze interval by 2 minutes at a time