

External microphone

BE9127



Features

Use it during meetings
To improve speech intelligibility

Use it for TV-listening
For clear and distinct sound

5 meters long cable
No extension cables required


For meetings and TV-listening.

The external microphone is specially developed for the Maxi and Mino digital speech amplifiers. It can be used as a remote microphone during meetings or conferences to improve the speech intelligibility in busy environments.

The external microphone can also be used for listening to the TV. Just attach the microphone close to the TV's loudspeaker and enjoy the loud and clear TV-sound.

Technical specifications

| | |
|------------------|---|
| ■ Cable length | 4.9 m, 16.3' |
| ■ Connectors | 2.5 mm mono jack plug |
| ■ Directivity | Omni |
| ■ Weight | 35 g, 1.2 oz. |
| ■ Colour | Black |
| ■ In the package | External microphone Clip w. magic tape sticker Microphone foam cap Getting started guide |



© 2016, Bellman & Symfon AB. All rights reserved.
Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB. This leaflet is subject to change without prior notice.

BE9127_001DAS002

bellman.com