# Microset<sup>™</sup>

**BE9129** 



### Bellman & Symfon<sup>®</sup>



**Features** 

Hands-free listening Tuck away Mino in your pocket

Small and lightweight For a discreet look

Works out of the box No fitting procedure required

## All-in-one microphone and earbud for the Mino.

50 Hz to 20 kHz

10 mW

32 ohm

The MicroSet is specially designed for the Mino speech amplifier. Put the Mino in your inner pocket and use the lapel microphone to pick up the sound. It includes an omnidirectional lapel microphone and a single earphone.

#### **Earphone specifications**

- 113(± 3) dB/mW (@ 1kHz) Sensitivity
- Frequency response
- Maximum power
- Impedance

#### **Microphone specifications**

- Electret condenser microphone Type
- Omnidirectional Directivity
- Sensitivity
- Frequency response
- S/N ratio
- Bias voltage Max. 10 VDC

#### **Cable and connector**

- Cable length
- Connector

95 cm, 3.1' 3.5 mm stereo jack plug

-47 (±3) dB1

>58 dB

110 Hz to 16 kHz



#### © 2016, Bellman & Symfon AB. All rights reserved.

Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB. This leaflet is subject to change without prior notice.

BE9129\_001DAS002

