

Maxi Personal amplifier



Discover the joy of better hearing

Maxi is a personal amplifier that will help you enjoy the company of your near and dear ones. It brings out the words so you can hear them loud and clear, even when there is a lot going on in the background. This means you don't have to make so much of an effort to hear.

Use Maxi for example when the family gathers around the dinner table, in the store, or to listen to the radio and TV. Maxi delivers clear and distinct speech and works both with and without hearing aids.



**”Now I am not
so tired after the
family dinner”**





” Now we no longer talk past each other ”

At home

Even in quiet surroundings, it can sometimes be difficult to distinguish the details of a conversation.

Just place Maxi on the table in front of the person you are talking to and it will take care of the rest. You can choose to listen with headphones, earphones or Stetoclips. Maxi is also an excellent complement to hearing aids when there is a lot going on in the background. It will filter out the background noise and clarify speech further before it reaches the hearing aid.

- ▼ Headphones
- ▼ Earphones
- ▼ Stetoclips



At social events

Family parties are often noisy events with high sound levels and a lot of people talking at the same time. On such occasions it can be particularly challenging to follow a conversation. After tiring of having to ask people to repeat themselves, it is easy to end up silent.

With Maxi you can hear the words even in busy environments. Just point it to the person talking and it will make the words stand out from the background noise.

- ▼ Neck loop

A neck loop is used to connect Maxi to a hearing aid or a cochlear implant.





By the TV

If you have difficulty hearing, it is easy to turn the volume on the TV and radio to the max. This is not always appreciated by the other family members.

Now there's no reason to turn up the volume to hear better. Connect Maxi to the TV or radio with the accompanying cable and you can enjoy great sound in stereo while sparing your family's ears. If you want to comment on the show, simply press the Push-to-Listen button to join the conversation.

▼ Cable kit, 16 ft (5 m)



In theaters

A visit to a movie or concert is a pleasant break from everyday life. But it is easy to miss out on the experience if you have to make an effort to hear.

Many public venues, such as theaters, concert halls, and churches, are fitted with loop systems. In these venues, Maxi captures the sound from the orchestra, movie or speaker with the built-in Telecoil. With a simple press of a button, you can switch to the internal microphone when you want to talk to the person next to you.

▼ T-symbol

Look for a sign with a T-symbol next time you visit your favorite movie theater or concert hall.





”Now I don’t need to worry about getting it wrong”

At the pharmacy

For people who are hard of hearing, something as simple as a visit to the pharmacy can become a real challenge. When the pharmacist’s voice blends in with the background it is easy to miss out on important information. On the other hand, you don’t want your private conversation to be overheard by others.

Place Maxi on the counter to hear the advisor loud and clear, without her having to raise the voice. You can also bring it to your next doctor’s appointment.

In the car

A car is one of the toughest environments to conduct a conversation in. Besides the noise from the traffic, the passengers are normally facing away from you. This can make it very hard to hear what’s being said.

Just hang Maxi around the neck of your travel companion and it will make the words stand out from the traffic noise. If you are sitting further away from each other you can connect a lapel microphone that you attach to your travel companion’s sleeve with a clip.

▼ Travel case



Do you travel occasionally? Order the travel bag in durable canvas, specially designed to protect Maxi and the accessories.

▼ External microphone



Easy to use

Up and running within minutes

Maxi is intuitively easy to use and you are up and running within minutes flat. Everything is easily accessible with the simple press of a button and there are no complicated menus to learn. Just connect your earphones or neck loop and adjust the volume to your liking



Great ergonomics


High quality and robust

Maxi is very robust and user friendly. It features large tactile buttons, soft grip materials and high contrast indicators that make it very easy to handle. Maxi is a great choice for those with low vision or limited dexterity.

Clear and distinct sound

Listen without effort

Maxi is equipped with the latest digital technology and delivers clear and distinct speech of the same quality as modern CD and DVD players. It efficiently amplifies speech while filtering out unwanted background noise. This means, you can take part in the conversation without having to put so much effort into just getting the words right.



Using Maxi with hearing aids

Maxi clarifies speech and filters out disturbing background noise before it reaches the hearing aid or cochlear implant. This means you don't have to make so much of an effort to get the words right. Here is how it works.



1. On the table

Maxi captures the sound, filters out the background noise and clarifies the speech.

2. Around the neck

The neck loop transmits the speech wirelessly to the hearing aid or cochlear implant.

3. In the ear

The hearing aid or cochlear implant refines the sound further before it reaches the ear.

Requires a hearing aid with Telecoil



Features

- Omnidirectional microphone ✓
- Digital signal processing ✓
- Digital noise reduction ✓
- Tone control ✓
- Volume control ✓
- Balance control ✓
- Built-in Telecoil ✓
- Input for external sound source ✓
- Push-to-Listen ✓
- Settings menu ✓

In the box

- Maxi personal amplifier
- Neck strap
- Belt clip
- User manual
- 2 x AA alkaline batteries

Content may vary with territory



Dedicated to your hearing

Bellman & Symfon of Sweden is dedicated to improve the quality of life for people with hearing problems. Our people and partners are devoted to this mission and we work closely with hearing healthcare professionals and leading experts to protect your hearing and make better hearing possible for everyone.

European office

Bellman & Symfon Europe AB
Södra Långebergsgatan 30
421 32 Västra Frölunda
Sweden
Phone +46 31 68 28 20
Fax +46 31 68 28 90
E-mail info@bellman.se

Asia Pacific office

Bellman & Symfon Asia Ltd.
1115-16 , Main Tower, GD Int'l Hotel
339 Huanshi Dong Rd
Guangzhou, China, 510095
Phone +86 20 8318 2950
Fax +86 20 8333 5528
E-mail info.cn@bellman.com