

# Earphones

BE9124-P02



### Features

**Closed drivers**  
Minimizes the risk of feedback

**In-ear design**  
Increased speech intelligibility

**Small and lightweight**  
For a discrete look

## Lightweight and discreet Earphones.

The Earphones are specially selected for the Bellman & Symfon range of digital speech amplifiers and listening systems. They offer high-clarity stereo sound with great detail and feature closed drivers that minimize the sound leakage and risk of feedback.

The Earphones are lightweight and provide excellent isolation from ambient noise, increasing the overall speech intelligibility.

### Technical specifications

- |                                  |                      |
|----------------------------------|----------------------|
| ■ Frequency range                | 20 Hz - 20 kHz       |
| ■ Impedance                      | 16 $\Omega$ +/- 10%  |
| ■ Sensitivity                    | 93 dB/mW @ 1 kHz SPL |
| ■ Connector                      | 3.5 mm stereo plug   |
| ■ Cable length                   | 1.25 m, 49"          |
| ■ Weight                         | 16 g, 0.6 oz         |
| ■ Colour                         | Black                |
| ■ Tangle-free cable              |                      |
| ■ Two sizes of ear tips included |                      |

### Approvals

**MD** This product is a medical device in conformity with Regulation (EU) 2017/745 (Medical Device Regulation).

© 2025 Bellman & Symfon AB. All rights reserved.  
Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB. This leaflet is subject to change without prior notice.

BE9124-P02\_001DSH1.0-EN [bellman.com](https://bellman.com)