

Technical specifications

Neckloop diameter: 22 cm, 9"

Weight

BE9159: 62 g, 2.2 oz BE9161: 58 g, 2.0 oz

Cable length

BE9159: 90 cm, 3' BE9161: 15 cm, 6"

Connectors

3.5 mm tele plug (stereo) gold plated connector, 90 degrees angle
(Break-away connector on cable)

Load impedance: 2 x 5 Ω

Magnetic output

1500 mA/m @ 15 cm, 6" distance and
2 x 50 mW input signal

Operation, storage and transport conditions

Operating temperature: 0°–35° C, 32°–95° F.

Transport and storage temperature: -10°–50° C,
14°–122° F.

Relative humidity: 15% – 90% (non condensing).

North America service

Bellman & Symfon North America Inc.
5509 Business Dr, Unit B
Wilmington, NC 28405
Unites States
Phone: +1 (910) 239 5796
Toll-free: + (833) 235 5626
E-mail: support.bsn@bellman.com

Asia Pacific service

Bellman & Symfon Asia Ltd.
1401-06, Main Tower, GD Int'l Hotel
339 Huanshi Dong Rd
Guangzhou, China, 510095
Phone +86 20 8318 2950
Fax +86 20 8333 5528
E-mail info.cn@bellman.com




Manufacturer

Bellman & Symfon Group AB
Södra Långebergsgatan 30
436 32 Askim Sweden
Phone +46 31 68 28 20
E-mail info@bellman.com
bellman.com



Revision: BE9159_006MAN010
Date of issue: 2020-09-15
TM and © 2020 Bellman & Symfon AB.
All rights reserved.




User information

- This Neckloop from Bellman & Symfon is to be used with hearing aids equipped with Telecoil.
- It is optimized for the Maxi, Sonido, Mino speech amplifiers and the Domino listening system.
- The neckloop can also be used with smart devices, TVs, radios or other equipment that features a 3.5 mm stereo headphone jack intended for low impedance headphones.
- It is discreet and comfortable, and can easily be worn under your shirt. It also incorporates a breakaway neck strap for increased safety.
- **Important!** If you have a pacemaker, we recommend you check with your general practitioner or cardiologist before using a neckloop.



Neckloop

BE9159/BE9161



Use with Domino or Mino

