

Listen at the volume of your choice

Connect Domino to the TV to enjoy great sound in stereo while sparing your family's ears. You can also listen to your HiFi system or follow your favorite show on the radio without disturbing a living soul. If you want to talk to the person next to you, just press the microphone button, so you can talk inhibited.



The Domino system consists of a transmitter and a receiver. This is how you use it.

The receiver captures speech close to you. Wear it around your neck or place it on the table in front of you.

The transmitter picks up speech from further away. Just place it across the table – or ask your colleague, schoolmate or teacher to wear it.



Dedicated to your hearing

Bellman & Symfon of Sweden is dedicated to improve the quality of life for people with hearing loss. Our people and partners are devoted to this mission and we work closely with hearing healthcare professionals and leading experts to protect your hearing and make better hearing possible for everyone.

Domino Pro Personal listening system



Discover the joy of better hearing

Domino is a listening system that helps people with hearing loss to hear better in busy situations, where hearing aids alone may not be enough. It brings out the words so you can hear them loud and clear, even when there is a lot going on in the background. This means you don't have to make so much of an effort to hear what people are saying.

Use it, for example, in a business meeting, when you're out dining with friends, or to hear a distant speaker at a seminar.



Hear without effort

Domino uses digital technology to enhance speech and filter out unwanted background noise before it reaches your hearing aid or speech processor. The clean and crisp sound makes it easier to hear the words so you can join the conversation without effort.

Just point and listen

The Domino units are equipped with two different microphones. While the omnidirectional microphone is useful in quiet surroundings, you can use the directional microphone to 'point' to the person talking in situations when it is particularly difficult to hear.

This is especially useful in busy environments, such as canteens, restaurants and terminals, where it would otherwise be difficult to distinguish speech.



European office

Bellman & Symfon Europe AB
Södra Långebergsgatan 30
421 32 Västra Frölunda
Sweden
Phone +46 31 68 28 20
Fax +46 31 68 28 90
E-mail info@bellman.se

Asia Pacific office







Bellman & Symfon Asia Ltd.
1115-16, Main Tower, GD Int'l Hotel
339 Huanshi Dong Rd
Guangzhou, China, 510095
Phone +86 20 8318 2950
Fax +86 20 8333 5528
E-mail info.cn@bellman.com

© 2013, Bellman & Symfon AB, all rights reserved. Bellman, Bellman & Symfon and the Bellman & Symfon logo are registered trademarks of Bellman & Symfon AB.

bellman.se/domino

BE2210_0928BH001_EN

 **Bellman & Symfon®**
DESIGN FOR EARS



How Domino works



- 1. On the table**
The receiver captures the sound nearby, while the transmitter picks up the sound further away. The sound is then streamed to the receiver that enhances speech and reduces disturbing background noise.
- 2. Around the neck**
The neck loop transmits the speech wirelessly to the hearing aid or cochlear implant's speech processor.
- 3. In the ear**
The hearing aid or speech processor refines the sound further before it reaches the ear.



* If you don't use hearing aids, you can choose to listen with headphones or earphones.

In the box

Everything included. The package comes complete and ready to use. Discover the joy of better hearing with the Domino Pro personal listening system.

- ✓ Domino Pro transmitter
- ✓ Domino Pro receiver
- ✓ Durable travel case
- ✓ In-ear earphones
- ✓ Tie-clip microphone
- ✓ Stereo cable, 1.5 m
- ✓ Neck straps
- ✓ Universal charger that works on all continents
- ✓ User manual
- ✓ Pocket reference guide

